




SNAP OFF SCREW SYSTEM

Self-drilling and Self-tapping

HBS COMPRESSION SCREW SYSTEM

Double threading for optimal compression

Scaphoid Bone Screws TITAN



HAND AND FOOT SURGERY



INSTRUMENTS



Recommended Sterilization Container for Set
Cat. No. T4700 (are not include in the Set)
Cat. No. 15-2162 Sterilization Container 310x190x65 mm,
Lid and Bottom perforated



Cat. No. T4701
Graphic Case (Tray in Steel) for Set Cat. No. T4701

INSTRUMENT / IMPLANT SET

Cat. No. T4700 HBS and Snap off Screw
Implant and Instrument Set



IMPLANTS

For Cat. No. see page 6!

Scaphoid bone screws in Titanium DIN ISO 5832-3
HBS bone screw system

Dia. 4.0 mm  Dia. 3.0 mm

For Cat. No. see page 9!

Self-drilling and tapping Snap off Screw Dia. 2.0 mm
in Titanium DIN ISO 5832-3 (weil osteotomy)

 Dia. 2.0 mm

Listing for Set Cat. No. T4700

| Cat. No.: Instruments for HBS - Snap off Screw System | Pcs.: |
|---|-------|
| 4-739 HBS - Guide Wire Dia. 1.0 mm x Length 80 mm | -4- |
| 4-740 HBS - Cannulated Drill Bit Dia. 2.1 mm / 3.3 mm | -1- |
| 4-741 HBS - Cannulated Drill Bit Dia. 2.1 mm / 3.3 mm | -1- |
| 4-742 HBS - Screw Length Gauge | -1- |
| 4-743 HBS - Measuring Sleeve for Guide Wire | -1- |
| 4-744 HBS - Cannulated Screw Driver, Hexagonal 2.0 mm | -1- |
| 4-745 HBS - Screwdriver for Snap off Screws | -1- |
| 4-184 Screw Forceps | -1- |
| Cat. No.: Dia. 3.0 mm Ti. HBS Screws, Cannulated | Pcs.: |
| 4-712 Length 12 mm, Titanium | -4- |
| 4-714 Length 14 mm, Titanium | -4- |
| 4-716 Length 16 mm, Titanium | -4- |
| 4-718 Length 18 mm, Titanium | -4- |
| 4-720 Length 20 mm, Titanium | -4- |
| 4-722 Length 22 mm, Titanium | -4- |
| 4-724 Length 24 mm, Titanium | -4- |
| 4-726 Length 26 mm, Titanium | -4- |
| 4-728 Length 28 mm, Titanium | -4- |
| 4-730 Length 30 mm, Titanium | -4- |
| Cat. No.: Dia. 2.0 mm Ti. Snap off Screws | Pcs.: |
| 4-760 Length 11 mm, Titanium | -4- |
| 4-762 Length 12 mm, Titanium | -4- |
| 4-764 Length 13 mm, Titanium | -4- |
| 4-766 Length 14 mm, Titanium | -4- |



INSTRUMENTS



Cat. No. 4-739 HBS- Guide Wire Dia. 1.0 mm x Length 80 mm



Cat. No. 4-740 HBS- Cannulated Drill Bit Dia. 2.1 mm/3.3 mm



Cat. No. 4-741 HBS- Cannulated Drill Bit Dia. 2.1 mm



Cat. No. 4-742 HBS- Screw Length Gauge



Cat. No. 4-743 HBS- Measuring Sleeve for Guide Wire 4-739



Cat. No. 4-744 HBS- Cannulated Screw Driver, Hex. 2.0 mm



Cat. No. 4-745 HBS- Screw Driver for Snap off Screws



Cat. No. 4-184 Screw Forceps



9-540 Verbrügge Bone Holding Forceps 175 mm



Cat. No. 13-3250 Inge 16.0 cm




DOUBLE THREADING FOR OPTIMAL COMPRESSION

Design rationale & Main features

The Compression Screw is easy to insert (over a guide wire) and provides efficient compression (through two separate threadings with different pitches, and an intermediate unthreaded section), thus ensuring quick, dependable internal fixation.

INDICATIONS

- Distal & proximal metatarsal osteotomies
- SCARF osteotomy
- Uni and biocortical internal fixation (ex.: scaphoid)
- Small bone fusion



| Cat. No.: | Total Length | Cat. No.: | Total Length |
|-----------|--------------|-----------|----------------|
| | Cannulated | | non Cannulated |
| 4-712 | 12 mm | 4-713 | 12 mm |
| 4-714 | 14 mm | 4-715 | 14 mm |
| 4-716 | 16 mm | 4-717 | 16 mm |
| 4-718 | 18 mm | 4-719 | 18 mm |
| 4-720 | 20 mm | 4-721 | 20 mm |
| 4-722 | 22 mm | 4-723 | 22 mm |
| 4-724 | 24 mm | 4-725 | 24 mm |
| 4-726 | 26 mm | 4-727 | 26 mm |
| 4-728 | 28 mm | 4-729 | 28 mm |
| 4-730 | 30 mm | 4-731 | 30 mm |



SELF-DRILLING & SELF-TAPPING

Design rationale & Main features


The Snap off Screw provides superior fixation: it saves time (no need for a pilot drill hole), and it is easy to use, safe (clean break), accurate (guide point), and efficient (self-compression).

The Snap off Screw consists of two parts: implantable screw which provides firm anchorage inserter which allows powered insertion.

INDICATIONS

- Weil osteotomy
- Unicortical internal fixation

- 4-760 Snap off Screw Ø 2 mm Length 11 mm
- 4-762 Snap off Screw Ø 2 mm Length 12 mm
- 4-764 Snap off Screw Ø 2 mm Length 13 mm
- 4-766 Snap off Screw Ø 2 mm Length 14 mm



SURGICAL TECHNIQUE (WEIL OSTEOTOMY)

EXPOSURE

The procedure is performed using a dorsal intermetatarsal and/or transverse approach.

After the two extensor muscles have been separated:

- Hohmann retractors are placed on both metatarsal sides.
- The metatarsophalangeal joint is dislocated between the extensor digitorum longus and the extensor digitorum brevis.
- A Hinge spreader is inserted to protect the extensor muscles and afford good exposure for the osteotomy.





TREU = INSTRUMENTE

OSTEOTOMY

Osteotomy is performed using an oscillating saw:

- Make a 3 cm (approximately) horizontal cut parallel to the sole, to increase the interfragmental contact area and thus enhance healing.
- Osteotomy results in spontaneous recession of the metatarsal head, which relieves tension on soft tissue.

2



TRANSLATION

- Grasp the metatarsal head with Kocher forceps
- Use the „Index Plus Minus“ formula and the Lelièvre Curve to determine the amount of recession of the metatarsal head.
- The metatarsal head must be held in the correct position for subsequent screw fixation.

3



INSERTION OF THE SNAP OFF SCREW

- Connect the screw inserter to the power drill, and drive the screw into the metatarsal.
- The inserter snaps off as soon as the screw head makes contact with the dorsal cortex.
- If necessary, insertion of the screw can be completed with the special screwdriver (with 3 notches).


4



RESECTION OF THE BONE PEAK

- Bone peak is resected using Liston pliers. This allows deep flexion of the metatarsophalangeal joint.
- It may be necessary to perform a Z-shaped release (Green technique) of the extensor muscles.

5



© TREU GmbH / Atelier Urban A1 / 07/05

TREU Instrumente GmbH · Altentalstraße 6 - 10
D - 78532 Tuttlingen · Germany



Tel. +49 (0) 7461 / 9621-0 · Fax +49 (0) 7461 / 77567
info@treu.com · www.treu.com