- D** HANDHABUNGSHINWEIS
- GB** INSTRUCTIONS FOR USE
- F** INSTRUCTIONS DE SERVICE
- I** ISTRUZIONI PER L'USO

## Elektrodenhandgriff

## Electrosurgical Pencil

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100-001-040 bis 100-279-316

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**English Electrosurgical Pencil**

**1 Attention**

Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments.

**2 Intended Use**

The product must be connected with the appropriate output of the electrosurgical generator and serve to connect the generator with monopolar electrodes.  
Only to be used by skilled medical professional who have been introduced to the product, f.e. by this IFU.

Maximum frequency fmax: 4,0MHz  
Maximum rated voltage of accessory Umax: 4,3 kVp

In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used. (See also section "Intended Use").  
See this Instructions for Use, or the current product catalog for the maximum rated voltage of the product.  
In case of uncertainties, contact the manufacturer.  
Read the entire IFU of this product and of any used accessory as well as HF-generator and HF-neutral electrode (monopolar application).  
All requirements, safety notices and warnings included in the respective IFU's have to be followed strictly.  
This medium is not intended for users in the USA.

**Indication:**  
Product intended for open or endoscopic surgery and is used as contact and switching element.

**Contraindication:**  
Product is not intended for direct contact with the heart or the central circulatory or central nervous system.  
Application of high-frequency current may interfere with cardiac pacemakers and in vivo heart defibrillators, so affected patients must consult a cardiologist prior to the intervention.  
It is not allowed the patient contacts any grounded metal objects during the intervention.

**Connection and Activation:**  
Prior to application of monopolar products, it has to be safeguarded that an HF-neutral electrode has been applied correctly to the patient and that this device is connected properly to the respective HF-generator.  
Connect to the intended output of the HF-generator by the means of an appropriate cable.

Follow instructions set forth in the User Manual of the HF-generator and follow the aforementioned defined performance limitations (Umax and frequency).  
**Activation is done by the means of the footswitch (HF-Generator) or by the means of the buttons of the product.**  
Blue button: COAG / Yellow button: CUT  
Only connect product to a switched-off HF-generator or to an HF-generator in Standby-Mode (if applicable)

**Combination / Compatibility:**

Prior to use, compatibility of the product with the intended HF-accessory and HF-generator has to be verified.  
In case of uncertainties, contact the manufacturer of this product or the manufacturer of the used accessory or HF-generator.  
Only use in combination with compatible products and connectors.  
Never connect to different or unknown power sources!  
The instructions, safety and warning notes in the instructions for use of the accessories and the HF generator used must be followed.

**Feasible HF-generators for aforementioned article:**  
According to the selected HF cable (plug).

■ ERBE ICC / ACC; T-Series	■ Karl Storz
■ Martin	■ Berchold
■ Integra	■ Ellman
■ Valleylab	■ EMED / ARTro200

Follow the instructions, safety notices and warnings included in the User Manuals or IFU's of the used accessories and used HF-generator.

**IMPORTANT:**  
Handle with utmost care.  
This does not only apply for the duration of the surgery but also for the complete duration of storage, processing and transport as well as during the process of connecting the product with the HF-accessory and HF-generator.  
This applies especially for the thin components and other sensitive areas, f.e. the insulation.

**Prior to first use** as well as **after each use**, the product has to be cleaned, disinfected and sterilized (see section "(Re-) Processing: Cleaning, Disinfection and Sterilization"). Improper use immediately will result in loss of warranty.  
Liability for any damages incurred will not be accepted.

**3 Safety Notices - WARNING**

See these instructions for use, the label, or the current product catalog for the maximum rated voltage of the product.  
Follow the aforementioned allowed frequency of HF current.

Prior to each application, the product has to be cleaned, disinfected and sterilized according to a validation procedure (EN ISO 17665).  
Prior to usage a visual inspection and function test has to be done (see section "Visual Inspection and Function Test").  
In case damages, deformation or similar is detected on the product, it is not allowed to use the device.  
It has to be replaced by a new one.

At least one (1) cleaned, disinfected and sterilized backup product has to be available.  
It is the responsibility of the user to determine the appropriate product size and product type according to their professional judgement and based on the patient's specific indication, preferred surgical technique and history, etc.

Prior to application, verify that the product is connected properly and firmly to the HF-generator and that the HF-accessory used also is inserted firmly and correctly.  
Regarding HF-Handles with smoke evacuation, makes sure the cap and the smoke evacuation tube area firmly attached and the smoke evacuation tube is connected properly to the smoke evacuation system.  
This must be done carefully, in order to avoid damages on the product and/or injuries of patient, surgical personnel or third party.  
Excessive force can damage the product.  
Therefore, the product has to be observed during the complete application.  
Damp insulate connectors may cause flashovers that could damage the connector or its insulation.  
Prevent fluids from penetrating the connector during use.

**Exclusion:**  
To plug in and unplug product, grasp it only on the plug part.  
Do not activate the HF-generator during connection of the product.  
Do not link the product.  
Do not use wound up or rolled up products.  
Do not activate the product as long as it is in contact with metal objects and/or optics.  
During an electrosurgical intervention the patient must not come into contact with grounded metal objects such as f.e. surgical desk frames, instrument trays etc.  
Pay attention that no flammable substances are present in immediate vicinity as otherwise there is danger of explosion.

**(Refer to section "General Safety Notices for HF-Technology")**  
After switching off the electrosurgical current, the product may still be hot and can lead to unintentional burns.  
Do not activate the product uninterruptedly during a longer period.  
It is the responsibility of the user to apply low performance settings on the HF-Generator in order to achieve the envisaged effect for the respective intervention.

During the surgery, mechanical forces may lead to deformation or signs of wear of the product.  
The application of electrostimulatory current can lead to damage of heart pacemakers and in vivo heart defibrillators.  
Therefore, so affected persons have to consult a cardiologist prior to the intervention.

**4 General Safety Notices for HF-Technology (Excerpt)**

In addition to the acknowledged benefits of HF-Technology, the application includes several risks that have to be attended:  
Improper use and non-observation of the IFU can lead to unintentional burns of the patient, user or third party.  
Continued further education of the surgical personnel is recommended.

**a) Environment**  
Pay attention that no flammable substances (anesthetics, oxidizing gases, endogenous gases etc.) are present in immediate vicinity as otherwise there is danger of explosion.  
Only use non-flammable cleaning or disinfection agents.  
All oxygen connections must be tight and leak-proof during the procedure.

**b) Patient Positioning and Preparation**  
Ensure proper patient positioning, i.e., only use insulating surgical drapes that are dry, absorbent and liquid-tight.  
Isolate conductive surfaces and contact points towards the patient.  
Dry tissues of pulp are necessary to use for skin and breast folds as well as between the extremities.

Prior to application, remove any liquids that potentially accumulated in body cavities.  
Only use non-flammable disinfectants.  
Do not use alcohol-based tinctures f.e. Only use non-conductive irrigation fluid, if medically possible.

Attend the requirements on irrigation fluid for monopolar and bipolar products.  
Prior to application, remove any body jewelry from the patient.  
Putting a band aid over the body jewelry is not sufficient!

**c) Connections**  
Prior to application, make sure the product is connected properly to the HF-generator and make sure the correct power and performance setting is adjusted and displayed.  
It is the responsibility of the user to apply low performance settings on the HF-Generator in order to achieve the envisaged effect for the respective intervention.

**d) HF-Neutral Electrode for Monopolar Application**  
In case of monopolar application, select a HF-neutral electrode suitable for the patient, apply it correctly and connect it properly with the respective HF-generator.  
Follow all instructions for proper application of the HF-neutral electrode, incl. patient protection and patient monitoring, monitoring of the HF-neutral electrode and all further provisions, safety notices and warnings included in the IFU of the HF-neutral electrode.

**e) Patient Reactions**  
All electrosurgical devices potentially can cause muscle stimulation during the application.  
The design of this product minimizes the risk of this undesirable effect.  
Nevertheless, muscle stimulation can lead to an unexpected movement of the patient in the surgical field.

**f) Handling HF-Accessory**  
Make sure the accessory is compatible.  
Do not touch the instrument tip during the complete application.  
As long as the product is not applied, place it on a dry, clean and non-conductive and well-visible surface, that is not in contact with the patient.  
Never store product on the patient.  
Unintended activation of the product can lead to unintentional burns or other injuries of patient, user or third party.  
Never wrap cables around the patient and never lay cables over the patient.  
Apply only short activation times.  
Keep on longer breaks between activation phases.  
Only adjust low power settings.

**g) Completeness of the System**  
At the end of the surgery, confirm completeness of the system.

**5 (Re) Processing: Cleaning, Disinfection and Sterilization**

**5.1. Maximum number of reprocessing cycles**  
**Due to the design, materials used, intended use as well as wear and tear, a maximum limit of performable reprocessing cycles cannot be determined.**  
When applied according to the Intended Use, the product underlies natural wear and tear, considering manner and duration of the application as well as manner and frequency of reprocessing.  
Therefore, a visual inspection and function test has to be done prior to each usage. (Refer to section "Visual Inspection" and "Function Test").

Visual inspection and function test, especially the condition of the complete cable, insulation, port, plug component and pins are decisive for whether the product is allowed to be applied again.

**5.2. Time requirements for cleaning and disinfection.**  
**Preparation for cleaning, pre-cleaning and automated cleaning and disinfection has to be done immediately after the application, however not later than 1 hour after the application.**  
**Avoid idle time.**

**5.3. After the Application**  
Clean and disinfect product immediately after the application.  
However not later than 1 h after application.

After application of the product, deposit the product carefully (protection of lifetime of the product).  
After application, separate contaminated product and deposit it in a suitable container (deposit means "do not drop").  
Immediately remove gross stain.  
Immediately mark damaged or defective products.  
Accessory, that do not fit for the steves of the cleaning and disinfection device (CDD), shall be deposit separately in suitable containers.

Close firmly all disposal/transport containers, in order to avoid drying of stain.  
Organize transport of contaminated products the way that contamination of the transport ways and environment also is avoided (closed transport).  
Unused reusable products have to be reprocessed as well. (Refer to AKI (German organization for instrument reprocessing), Red Brochure, page 30-32).  
Take care, all and any transport containers are cleaned and disinfected after transport as well.

**5.4. Validation of (Re)Processing**  
The following validated processing procedure is recommended.  
Equivalent deviant processes are possible.  
Then, it is the sole responsibility of the user to safeguard the suitability of the actually applied procedure by suitable means (f.e. validation, routine monitoring, verification of material compatibility etc.).  
Automated cleaning and disinfection always are preferable.  
The following procedure was validated according to EN ISO 17665.

Additional applicable processing requirements specific to the respective clinical place (operator) as well as national or country specific regulations have to be followed as well.  
Never use sharp objects for cleaning.  
Disinfectants always have to be rinsed and removed carefully.

**5.5. Preparation for Cleaning**  
New product. Remove product from their packaging.  
Used product. Any attached HF-accessory (f.e. HF-Adapter, HF-Electrode etc.) has to be disconnected.  
Place the product in a container intended for pre-cleaning.

**5.6. Pre-Cleaning**  
The product has to be cleaned immediately after usage (see section above).  
Prepare a cleaning bath using 8 ml cleaning agent (f.e. CIDENZYM®) per 1 l water.  
Attend requirements set forth by the manufacturer of the cleaning agent.

**Requirements on cleaning agent:**  
Intended for medical devices made from metal and plastics.  
pH value: 5.5 to 12.3  
Do not use high alkaline cleaning agents.  
These will impair the lifetime of the product.  
Attend reaction time set forth by the manufacturer of the cleaning agent.  
Take care the product will not touch other parts in the cleaning bath.  
Thoroughly remove surface stain with a soft brush or synthetic fleece, as otherwise particles or dried secretions may adhere.

This could make subsequent cleaning and sterilization difficult or impossible.  
Do not use sharp, sharp edged or scrubbing devices for cleaning.  
Ensure that areas difficult to reach are cleaned thoroughly and rinsed several times.  
Rinse product thoroughly for at least 1 minute under running tap water, in order to remove residues of the cleaning agent.

Check product visually for remaining stain.  
In case stain is still visible, repeat the process of pre-cleaning as often and long as it needs until all stain is removed and stain is no longer visible.  
This pre-cleaning step always has to be done prior to the manual cleaning or cleaning with the cleaning and disinfection device (CDD).

**5.7. Manual Cleaning and Disinfection**  
Manual cleaning and disinfection is not recommended.  
**5.8. Automated Cleaning and Disinfection**  
Only use cleaning and disinfection devices (CDD) with proven efficiency according to EN ISO 15883.

Follow data of the manufacturer of the cleaning and disinfectant machine.

**Requirements on cleaning agent:**  
Intended for medical devices made from metal and plastics.  
pH value: 5.5 to 12.3  
Do not use high alkaline cleaning agents.  
These will impair the lifetime of the product.  
Recommendation: neodisher® mediclean forte (Dr. Weigert GmbH & Co. KG).

**Requirements on CDD:**  
Apply program for thermal disinfection.  
Follow instructions and data regarding program course and machine.  
Product has to be stored safely and protected against mechanical damages during automated cleaning and disinfection.

Do not clean together with sharp edged or pointed objects.  
Deposit product in a suitable rinsing basket.  
Take care cable is not kinked or squeezed.  
Follow data for loading of the cleaning and disinfection device (CDD).

**Cleaning and Disinfection Program**  
Start the program course with following features:  
Thermal disinfection  
▪ 5 to 10 minutes at 90°C to 93°C (tolerance according to EN ISO 15883-1 or -2) A0 > 3000  
▪ Final rinsing with deionized water  
▪ Product drying

At the end of program course, remove product and inspect it on remaining stain.  
In case of residues, repeat automated cleaning and disinfection step as long as it takes until stain is no longer present.  
Dry cavities and not sufficiently dried areas with sterile compressed air <2 bar.

Immediately after removal of product and immediately after additional drying on a clean place, put product in a single use sterilization packing (double packing) from paper or foil or put product in a sterilization container.  
Respect requirements for sterilization packaging according to EN ISO 11607 and EN 868.

**5.9. Sterilization**  
Only products that have been cleaned and disinfected are allowed to be sterilized.  
Only apply steam sterilization in autoclave (fractioned pre-cleaning with sufficient product drying) for this product.

**Adjust sterilization parameters:**  
▪ Minimum 134°C and maximum 137°C in saturated steam.  
▪ Holding time at least 5 Min. until max. 20 Min.  
▪ Drying in vacuum for at least 10 Min.  
▪ Sterilizer (Class B) according to valid national standards and regulation (f.e. EN 13060 or EN 285).

**Example:** Sterilizer Class B, manufacturer: Tuttnauer.  
Respect data of the sterilizer's manufacturer regarding load, handling and drying times.  
**Exclusion:** Do not apply hot air, EO-gas, Radiation or Plasma for sterilization, or any other sterilization method for this product.

**IMPORTANT:**  
Prior to usage, let product cool to room temperature.  
It is the sole responsibility of the user to maintain the sterile condition of the product after the sterilization process.  
In case the aforementioned chemicals and machines for cleaning, disinfection or sterilization are not available, it is the responsibility of the user to validate the procedure actually applied.  
Also, if a sterilization method other than described above is applied, this deviating procedure has to be validated by the user accordingly.

**5.10. Limitation of Reprocessing**  
The product life time is depending on wear and tear, handling, application time, damages as well as frequency of reprocessing.  
Therefore, a visual inspection and function test has to be done prior to each usage. (Refer to section "Visual Inspection" and "Function Test").  
Only an undamaged product is allowed to be reused.

**6 Visual Inspection and Function Test**

Never lay this product on the patient or in their direct vicinity.  
A product exhibiting damages, pressure points or questionable condition is not allowed to be used and has to be replaced by a new one.  
As for products with buttons, check the correct function.  
Blue button: COAG / Yellow button: CUT  
If buttons do not work, use footswitch.  
Note: If the cable conductor is broken or the electrical continuity of the cable is interrupted in any other manner it can cause arcing in the patient return system or in the active circuit.

This can lead to the patient suffering burns or the outbreak of a fire.  
If damaged, do not use these products.  
Prior to each usage, the complete product, especially the insulation has to be inspected for pressure points and damages.

Prior to usage an electrical continuity test has to be done. In case the product does not pass the electrical continuity test, it is not allowed to use the product any longer and has to be replaced by a new one.

#### 7 Exclusion of Repair and Modification

Unauthorized repairs are strictly prohibited. Unauthorized modification is strictly prohibited and lead to an immediate loss of manufacturer's warranty. It is not allowed to bend product or to manipulate the product in any other manner.

Especially product having a hook as active part is never allowed to be bent.

This could lead to severe injuries of patient, user or third party.

#### 8 Packaging, storage and Transport

Store in a clean and dry environment. Storage only in protective cases with individual compartments or individually sealed in foil.

Handle with utmost care during transport, cleaning, disinfection, maintenance, sterilization and storage (see also section "after the application").

This applies especially for the thin components and their sensitive areas, f.e. the insulation.

Do not store or transport the product together with sharp edged or pointed objects.

Maintenance of the sterile condition after sterilization process is the sole responsibility of the user.

info@micromed.com

#### 9 Warranty

MICROMED Medizintechnik GmbH only provides inspected and faultless products. Our products satisfy the highest quality standards. Liability and warranty are excluded for all products that have been modified in any way, not applied according to their Intended Use, or that have been handled or applied improperly, or in case of any other deviation from the instructions set forth in this IFU.

Furthermore, MICROMED Medizintechnik GmbH denies any liability for any accidental, intentional damage or for a damage or loss arising out of handling or application of the product.

Additionally, all liability and warranty are omitted in case our product was repaired by a company that has not been authorized by us.

Unauthorized repairs are strictly prohibited.

#### 10 Return

Returned products only are accepted, if they are marked as "hygienically safe" or "not contaminated" and have been packed safely for the transport.

Use our Return Form for returns.

#### 11 Disposal

Disposal of the product, their packaging material as well as any accessory has to be done according to the applicable country specific requirements, regulations and laws.

Additionally, the applicable requirements of the respective clinical place(s) in regards to disposal of medical devices have to be followed as well.

#### 12 Regulatory Remark

Report serious incidents with the product to us. If you are a user in the European Union, also report incidents to the responsible authority in your Member State.

#### 13 About these Instructions for Use

Throughout the period of usage these IFU must be kept freely accessible for the user.

For a current revision of this IFU, please contact our customer service.

Changes reserved. This medium is not intended for users in the USA.

#### Explanation of symbols

	Batch Code
	Medical Device
	Catalog number
	Unique Device Identifier
	GS1 standard
	Caution!
	Consult Instructions for Use
	CE-marking
	Manufacturer
	Production date
	Keep dry
	Keep away from sunlight

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