DEUTSCH

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ENGLISH
Bipolar and Monopolar Electrodes / Adapter / Accessory kits

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 Please read all information contained in this insert attentively. Incorrect handling and care, as well as misuse, can lead to premature wear of surgical instruments or risks to patients and users.

Intended Use
 Bipolar electrodes for open surgery

The bipolar electrodes (REF 127-110-115 to 127-149-170 and 170-001-000 to 171-016-000) serve cutting and coagulation of biological tissue. The fully assembled instrument (if assembly is needed) has to be – with the appropriate cable - to monopolar or bipolar output of an HF generator. Only the defined parameters has to be used.

Maximum output voltage of the generator. U_{max}:
 For monopolare 4 kVp , (depending on the type, see catalogue)

For bipolare 250 Vp , 500 Vp (depending on the type, see catalogue)

Appropriate connecting cables for bipolar electrodes for open surgery:
 Micromed bipolar cable REF 130-1x0-x00.

Bipolar and monopolar electrodes
 The bipolar electrodes (REF 127-110-115 to 127-149-170 and 170-001-000 to 171-016-000) are used to coagulate biological tissue.

The monopolar electrodes (REF 120-400-001 to 124-315-110) are intended for general surgery, for coagulating or cutting biological tissue.

The following electrodes (REF 125-001-322 to 125-004-065) are intended for coagulating or cutting biological tissue. The electrodes are intended for use in minimally invasive and especially laparoscopic surgical procedures. The instrument has to be inserted through a trocar sleeve with the appropriate diameter. The monopolar hysterectomy loop (REF 125-00x-xxx) electrode is intended to be used for removal of the uterus in supracervical hysterectomies. The fully assembled instrument (if assembly is needed) has to be connected- with the appropriate cable - to monopolar or bipolar output of an HF generator. When indicated, monopolar or accordingly bipolar coagulation or cutting current can be selectively applied.

Maximum output voltage of the generator. U_{max}:
 for bipolar electrodes for MIS:
 500 Vp , 300 Vp , 250 Vp (depending on the type, see catalogue)

for monopolar electrodes for MIS:
 2kVp

Suitable cables for bipolar electrodes for MIS:
 Micromed bipolar cable REF 130-1x0-x00.

Suitable cables for monopolar electrodes for MIS:
 Micromed Monopolar Cable REF 102-0x4-x00.

 Instruments for electrosurgery must only be used by persons who have been specially trained or instructed in this.

Contraindications
 - Do not use the instrument if, in the opinion of the attending physician, the risks to the patient outweigh the benefits.

Incidents that have been reported in connection with the use of electrosurgical systems
 - Unintended activation with resulting tissue injury in the wrong location and/or damage to the equipment.
 - Fire in connection with surgical drapes and other inflammable materials.
 - Alternating current paths leading to burns on spots where the patient or user comes into contact with components without insulation.
 - Explosions caused by sparks in the proximity of inflammable gases.
 - Perforation of organs. Sudden severe bleedings.

Use and safety instructions
 Non-observance of these use and safety instructions may lead to injuries, malfunctions or other unexpected incidents.
 - When using electrosurgery in patients with pacemakers or other active implants, special requirements apply (e.g. low HF-current, patient monitoring). In any case, a cardiologist or appropriate medical specialist must be consulted.

- Before initial use and any further use, all instruments must be completely cleaned, disinfected and sterilised and their function must be checked.
 - It is very important to check every surgical instrument for visible damage and wear, such as cracks, breaks or insulation defects before each use. In particular areas such as blades, tips, notches, locking and blocking devices, as well as all movable parts, insulations and ceramic elements must be checked carefully.
 - Never use any damaged instruments.
 - Never use the instruments in the presence of flammable or explosive substances.
 - When temporarily not in use, the instrument must be placed electrically insulated from the patient.
 - Activate electrosurgical current only if the contact areas are in full view and have good contact with the tissue that needs to be treated. Do not touch any other metallic instruments, trocar sleeves, optics or similar objects during use.
 - Observe the use and safety instructions of the manufacturer of the high-frequency surgical device.
 Applies for monopolar mode of operation: Ensure correct application of the neutral electrode on the patient; otherwise, there is a danger of burns.

Assembly
 For assembly and disassembly of the instrument follow the pictogram, which is available upon request.

Special instructions for the hysterectomy loop electrode
 - Insert the instrument through the cannula with the loop pulled inside the shaft.
 - Thoroughly check correct placement of the loop around the uterus prior to activation of the electrical current; the entire tissue surrounded by the loop will be cut.
 - After correct placement, maintain tension of the loop to avoid uncontrolled slip of the loop. Maintain some safety distance to the ligature.
 - During electrosurgical cutting, pull the loop continuously inside the shaft. Cutting will start at the uninsulated part of the loop wire.
 - The loop wire is a replaceable, sterilisable product for single use that should be replaced following each procedure.

Reprocessing
 Due to the product design, the materials used and the intended purpose, it is not possible to define a limit with regard to the maximum possible number of reprocessing cycles. The serviceable life of the instruments is determined by their function as well as by a careful handling.
 Instruments for electrosurgery are by their nature subject to increased wear depending on the type and time of use.

Preparation and transport
 Immediately after each use, clean the instruments with a soft brush under cold tap water until all visible contamination is removed. Do not use fixation agents or hot water (>40°C). Storage and transport of the instruments to the reprocessing location must take place in a sealed container.
 Complex instruments must be taken apart for cleaning and disinfection in accordance with pictogram.

Machine reprocessing
Cleaning
 Place the instruments in a basket on the insert module or on the inserts of the MIS module and start the cleaning process.
 1. Pre-rinse with cold water for 1 min
 2. Discharge
 3. Pre-rinse with cold water for 3 min.
 4. Discharge
 5. Wash at 55°C with a 0.5% alkaline or at 45°C with an enzymatic cleaning agent for 5 min.
 6. Discharge
 7. Neutralise with warm tap water (>40°C) and a neutralising agent for 3 min.
 8. Discharge
 9. Rinse with warm tap water (>40°C) for 2 min.
 10. Discharge

Disinfection
 Machine-operated thermal disinfection must be carried out under observation of the national requirements regarding the A0 value (see ISO 15883).

Drying
 Dry the outside of the instruments by carrying out a drying cycle of the cleaning/disinfection machine.
 If necessary, manual drying may additionally be carried out using a lint-free cloth. Dry cavities by blowing with sterile compressed air.

Manual reprocessing

Ultrasonic pre-cleaning
 1. The instruments are placed in an ultrasonic bath with 0.5% enzymatic cleaning detergent and treated with ultrasound for 15 minutes at 40°C/104°F.
 2. Remove the instrument and rinse them completely with cold water to remove the cleaning detergent.

Cleaning
 Prepare a cleaning bath according to the manufacturer's instructions.
 1. Rinse products with cold tap water (<40°C) until all visible contamination has been removed. Remove adhering dirt by using a soft brush.
 2. Place products in the prepared cleaning bath so that they are completely submerged. Observe residence time according to the manufacturer's instructions.
 3. Clean the instrument in the bath manually using a soft brush. Brush all surfaces several times.
 4. **The following step only applies to channels and the insides of tubes:** Push the brush into and out of the tubes at least six times. Rinse the tubes with DI water. Repeat the procedure.
 5. Rinse the products thoroughly with DI water to remove the cleaning agents without residue.

Disinfection
 Prepare a disinfectant bath according to the instructions of the disinfectant manufacturer. Place the instruments in the disinfectant bath and observe the specified residence time. Rinse the products very thoroughly with DI water to remove the disinfectant without residue.


Drying
 Manual drying is carried out using a lint-free cloth and sterile compressed air, in particular for drying cavities and channels.

Functional test and packaging
 Perform visual inspection for cleanliness and integrity; if required, perform an assembly and functional test.
 If necessary, repeat reprocessing until the instrument is visually clean. Packaging must comply with the ISO 11607 and EN 868 standards for packaging for sterilised instruments.

Sterilisation
 Sterilisation of the products with fractional pre-vacuum procedure (in accordance with ISO 13060 / ISO 17665) under observation of the respective national requirements.
 - 3 pre-vacuum phases with a pressure of at least 60 mbar.
 - Heating up to a sterilisation temperature of at least 132°C and at most 137°C
 - Exposure time: at least 3 min.
 - Drying time: at least 10 min.
 If contamination with prions (CJD) is suspected, differing national guidelines are to be followed and longer holding times (i.e. 15 min.) may apply.

Storage
 Sterilised instruments must be stored in a dry, clean and dust-free environment. The applicable national guidelines must be followed.

Repairs
 Never attempt to perform repairs yourself. Service and repair work must only be performed by persons trained and qualified accordingly. If you have any question regarding these matters, contact either the manufacturer or your medico-technical department.

 Defective products must complete the entire reprocessing process before being returned for repair.

Information on the validation of the reconditioning
 The following testing instructions, materials and equipment have been used for validation:

Cleaning agents (for machine use):
 Neodisher FA by Dr. Weigert (alkaline)
 Endozyme by Ruhof (enzymatic)
Cleaning agents (manual cleaning):
 Cidezyme, Enzol Enzym detergent, Johnson&Johnson
Disinfectants (manual disinfection):
 Cidex OPA, Johnson&Johnson
Neutralising agent:
 Neodisher Z by Dr. Weigert
Cleaning and disinfection device:
 Miele Desinfector G 7735 CD
 Miele insert module E 327-06
 Miele MIS module E 450

For details, see report.
 SMP GmbH # 01707011901 (machine cleaning)
 MDS GmbH # 135196-10 (manual cleaning, sterilisation)
 Nelson Labs # 200432706-02 (sterilisation)
 MDS GmbH Testbericht 084183-10 (sterilisation)


If the chemicals and machines described above are not available, the user has to validate the used process accordingly.

Handling
 During transport, cleaning, care, sterilisation and storage, all surgical instruments should be handled with maximum care. This applies particularly to blades, fine tips and other sensitive areas.

Disposal
 Disposal must be carried out in accordance with the respective applicable local and national laws and regulations.

Warranty
 Micromed Medizintechnik GmbH exclusively supplies tested and faultless products to its customers. All products are designed and manufactured to comply with maximum quality requirements. We refuse any liability for products which have been modified as compared to the original product, misused or handled or used improperly.

Explanation of symbols

	Batch Code
	Medical Device
	Catalog number
	Unique Device Identifier
	GS1 standard
	Caution!
	Consult Instructions for Use
	CE-marking with notified body
	Manufacturer
	Production date
	Keep dry
	Keep away from sunlight