

MicroPen[®]

for RF-Units of
Micromed
ERBE
Martin
Berchtold
Aesculap
Valleylab



Cable

- in various lengths 3m, 4m, 5m
- with various types of plugs for different manufacturers' RF surgical units

Handle and cable separable

- flexible handle designed to fit into different RF-generators
- interchangeable cables

Electrode handle MicroPen



- aesthetic design
- optimum handling
- long life

Push button

- lift of stroke can be felt easily
- click control
- high security by durable push button
- push button completely covered with insulating silicone

Electrodes


- hexagonal fit for electrodes measuring 4 mm and 3/32 inches in diameter
- circular fitting for 1,6 mm microelectrodes
- electrodes can be exchanged easily



You want to know more?

To get further information
please visit www.micromed.biz


MICROMED®



100-104-016

100-201-040

100-201-024

100-101-040

100-203-040

MICROMED MICROPEN

12

12 MONTH GUARANTEE*

Optimum synthetic substance

- high voltage strength up to 7500 Volt
- great number of sterilization cycles
- highly resistant to detergents

Electrode handle MicroPen® • Highest quality with 12 month guarantee • Electrode handle MicroPen® from Micromed®

The **MicroPen®** is an innovative electrode handle that can live up to highest demands. Due to optimum technical production and the use of top quality materials we have created a handle with maximum longevity as well as optimum handling comforts.

Handle and cable separable

By separating handle and cable the components can be replaced separately in case of defects. This lowers costs and is environmentally friendly.

Highly convenient

Sensibly located pushbuttons and the easily felt pressure point (on/off click) as well as the rough outside surface make the MicroPen especially user-friendly.

Cleaning and Disinfection

The MicroPen can be sterilized more than 250 times.


Highest compatibility

The MicroPen is compatible to a wide range of different manufacturers' surgical units, as for instance: Aesculap, Berchtold, Integra, Erbe, Martin and Micromed.

*every MicroPen® comes with a 1 year guarantee on purchase.



www.micromed.biz



Micromed Medizintechnik GmbH
Esslinger Str. 37
D-78532 Tuttlingen

Tel +49 7462 923250
Fax +49 7462 923246
info@micromed.biz