





YOUR RELIABLE PARTNER IN TRAUMA CARE

WE FOCUS THE FUTURE

FIDES® HUMERAL NAIL SYSTEM







FIDES® HUMERAL NAIL SYSTEM

- antegrade and retrograde insertion
- cannulated and non-cannulated nail available
- 4 locking options distally
- radiolucent proximal targeting device
- optional adjustable distal targeting device
- user friendly instrumentation
- color coded locking screws


- Antegrade- und Retrograde Einbringung
- kanülierte- und nicht kanülierte Nägel verfügbar
- 4 distale Verriegelungsoptionen
- strahlendurchlässiges proximales Zielgerät
- optionales, verstellbares Distalzielgerät
- anwenderfreundliches Instrumentarium nach dem Minimalprinzip
- farbliche Kennzeichnung der Verriegelungsschrauben



 Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
1-3800	6.5 mm	190 mm	1
1-3802	6.5 mm	205 mm	1
1-3804	6.5 mm	220 mm	1
1-3806	6.5 mm	230 mm	1
1-3808	6.5 mm	240 mm	1
1-3810	6.5 mm	250 mm	1
1-3812	6.5 mm	260 mm	1
1-3814	6.5 mm	270 mm	1
1-3816	6.5 mm	280 mm	1
1-3818	6.5 mm	295 mm	1
1-3820	6.5 mm	310 mm	1
1-3822	6.5 mm	325 mm	1
1-3830	7.5 mm	190 mm	1
1-3832	7.5 mm	205 mm	1
1-3834	7.5 mm	220 mm	1
1-3836	7.5 mm	230 mm	1
1-3838	7.5 mm	240 mm	1
1-3840	7.5 mm	250 mm	1
1-3842	7.5 mm	260 mm	1
1-3844	7.5 mm	270 mm	1
1-3846	7.5 mm	280 mm	1
1-3848	7.5 mm	295 mm	1
1-3850	7.5 mm	310 mm	1
1-3852	7.5 mm	325 mm	1
1-3860	9.5 mm	190 mm	1
1-3862	9.5 mm	205 mm	1
1-3864	9.5 mm	220 mm	1
1-3866	9.5 mm	230 mm	1
1-3868	9.5 mm	240 mm	1
1-3870	9.5 mm	250 mm	1
1-3872	9.5 mm	260 mm	1
1-3874	9.5 mm	270 mm	1
1-3876	9.5 mm	280 mm	1
1-3878	9.5 mm	295 mm	1
1-3880	9.5 mm	310 mm	1
1-3882	9.5 mm	325 mm	1

FIDES® humeral nail (titanium), non-cannulated
FIDES® Humerusnagel (Titan), nicht kanüliert



 Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
1-3900	8.0 mm	190 mm	1
1-3902	8.0 mm	200 mm	1
1-3904	8.0 mm	210 mm	1
1-3906	8.0 mm	220 mm	1
1-3908	8.0 mm	230 mm	1
1-3910	8.0 mm	240 mm	1
1-3912	8.0 mm	250 mm	1
1-3914	8.0 mm	260 mm	1
1-3916	8.0 mm	270 mm	1
1-3918	8.0 mm	280 mm	1
1-3920	8.0 mm	290 mm	1
1-3922	8.0 mm	300 mm	1
1-3924	8.0 mm	310 mm	1
1-3926	8.0 mm	320 mm	1
1-3930	9.0 mm	190 mm	1
1-3932	9.0 mm	200 mm	1
1-3934	9.0 mm	210 mm	1
1-3936	9.0 mm	220 mm	1
1-3938	9.0 mm	230 mm	1
1-3940	9.0 mm	240 mm	1
1-3942	9.0 mm	250 mm	1
1-3944	9.0 mm	260 mm	1
1-3946	9.0 mm	270 mm	1
1-3948	9.0 mm	280 mm	1
1-3950	9.0 mm	290 mm	1
1-3952	9.0 mm	300 mm	1
1-3954	9.0 mm	310 mm	1
1-3956	9.0 mm	320 mm	1
1-3960	10.0 mm	190 mm	1
1-3962	10.0 mm	200 mm	1
1-3964	10.0 mm	210 mm	1
1-3966	10.0 mm	220 mm	1
1-3968	10.0 mm	230 mm	1
1-3970	10.0 mm	240 mm	1
1-3972	10.0 mm	250 mm	1
1-3974	10.0 mm	260 mm	1
1-3976	10.0 mm	270 mm	1
1-3978	10.0 mm	280 mm	1
1-3980	10.0 mm	290 mm	1
1-3982	10.0 mm	300 mm	1
1-3984	10.0 mm	310 mm	1
1-3986	10.0 mm	320 mm	1

FIDES® humeral nail (titanium), cannulated
FIDES® Humerusnagel (Titan), kanüliert





Cat. no. Art. Nr.			PU VPE
1-20250			1


FIDES® plug for humeral nail (titanium), cannulated
Verschlusskappe für FIDES® Humerusnagel (Titan)



Cat. no. Art. Nr.	Ø	Lenght Länge	PU VPE
1-23000	4 mm	10 mm	1
1-23020	4 mm	12 mm	1
1-23040	4 mm	14 mm	1
1-23060	4 mm	16 mm	1
1-23080	4 mm	18 mm	1
1-23100	4 mm	20 mm	1
1-23120	4 mm	22 mm	1
1-23140	4 mm	24 mm	1
1-23160	4 mm	26 mm	1
1-23180	4 mm	28 mm	1
1-23200	4 mm	30 mm	1
1-23220	4 mm	32 mm	1
1-23240	4 mm	35 mm	1
1-23260	4 mm	40 mm	1
1-23280	4 mm	45 mm	1
1-23300	4 mm	50 mm	1
1-23320	4 mm	55 mm	1
1-23340	4 mm	60 mm	1

FIDES® proximal locking screw 4.0mm (titanium)
FIDES® proximale Verriegelungsschraube 4.0mm (Titan)



	Cat. no. Art. Nr.	∅	Length Länge	PU VPE
	1-26000	3.5 mm	10 mm	1
	1-26020	3.5 mm	12 mm	1
	1-26040	3.5 mm	14 mm	1
	1-26060	3.5 mm	16 mm	1
	1-26080	3.5 mm	18 mm	1
	1-26100	3.5 mm	20 mm	1
	1-26120	3.5 mm	22 mm	1
	1-26140	3.5 mm	24 mm	1
	1-26160	3.5 mm	26 mm	1
	1-26180	3.5 mm	28 mm	1
	1-26200	3.5 mm	30 mm	1
	1-26220	3.5 mm	32 mm	1
	1-26240	3.5 mm	35 mm	1
	1-26260	3.5 mm	40 mm	1
	1-26280	3.5 mm	45 mm	1
	1-26300	3.5 mm	50 mm	1
	1-26320	3.5 mm	55 mm	1
	1-26340	3.5 mm	60 mm	1

FIDES® distal locking screw 3.5mm (titanium)
FIDES® distale Verriegelungsschraube 3.5mm (Titan)




	Cat. no. Art. Nr.		PU VPE
	1-20451	optional	1


adjustable distal targeting device for FIDES® humeral nail
verstellbares distales Zielgerät für FIDES® Humerusnagel



optional - not included in the set

 Cat. No. Art. Nr.	Set Cat. No. T 1002 Included Items	QTY St.
T 1006	tray for instrumentation set	1
1-20450	proximal targeting device (carbon)	1
1-20460	sleeve Dia. 6.0mm for 3.5/4.5mm screws	1
1-18640-0	calibrated drill bit dia. 2.7mm / 210mm	1
1-18640-2	calibrated drill bit dia. 3.2mm / 210mm	1
9-086	universal chuck with key cann. 4.1mm	1
1-2056	screw driver f. locking screws 3.5/4.0mm	1
1-1908	cannulated reamer Dia. 9.0mm	1
1-1876	tip threaded guide pin Dia. 3.2mm x 305	1
1-2040	nail guide wire Dia. 2.4 x 500mm	1
1-1882	open end wrench 17/14"	1
1-3718	driver for nail	1
1-37011	slide hammer for humeral nail, cannulated	1
1-1098	diamond point awl	1
1-20470	sleeve dia. 3.3mm for drill bit 3.2mm	1

instrument set for FIDES® humeral nail titanium
 Instrumentenset für FIDES® Humerusnagel



Cat. no. Art. Nr.			PU VPE
1-20450			1

proximal targeting device (carbon parts radiolucent)
Proximales Zielgerät (Carbonteile strahlendurchlässig)



Cat. no. Art. Nr.		Length Länge	PU VPE
1-20460		114mm	1
1-20460-1		130mm	1

sleeve dia. 6.0mm for 3.5/4.5mm screws
Hülse Ø 6.0mm für 3.5/4.5mm Schrauben



Cat. no. Art. Nr.		Length Länge	PU VPE
1-20470-1	distal	136mm	1

sleeve dia. 2.8mm for drill bit 2.7mm
Hülse Ø 2.8mm für Bohrer 2.7mm



Cat. no. Art. Nr.		Length Länge	PU VPE
1-20470	proximal	120mm	1

sleeve dia. 3.3mm for drill bit 3.2mm
Hülse Ø 3.3mm für Bohrer 3.2mm



Cat. no. Art. Nr.		Length Länge	PU VPE
1-18640-0	distal	210mm	1
1-18640-1	distal	230mm	1

calibrated drill bit dia. 2.7mm
kalibrierter Bohrer Ø 2.7mm



Cat. no. Art. Nr.		Length Länge	PU VPE
1-18640-2	proximal	210mm	1

calibrated drill bit dia. 3.2mm
kalibrierter Bohrer Ø 3.2mm




Cat. no. Art. Nr.			PU VPE
9-086			1

universal chuck with key cannulated. 4.1mm
Universalbohrfutter mit Schlüssel, kanüliert 4.1mm



Cat. no. Art. Nr.	hex sw		PU VPE
1-2056	2.5		1

screw driver for humeral locking screws 3.5/4.0mm
Schraubendreher für Humerus Verriegelungsschraube 3.5/4.0mm



Cat. no. Art. Nr.			PU VPE
1-1908			1

cannulated reamer dia. 9.0mm
kanülierter Fräser Ø 9.0mm



Cat. no. Art. Nr.			PU VPE
1-1876			1

tip threaded guide pin dia. 3.2mm x 305mm
Führungsdraht mit Gewinde Ø 3.2mm x 305mm



Cat. no. Art. Nr.			PU VPE
1-2040			1

nail guide wire dia. 2.4 x 500mm
Führungsdraht Ø 2.4 x 500mm für Nagel



Cat. no. Art. Nr.			PU VPE
1-1882			1

open end wrench 17/14
Gabelschlüssel 17/14



	Cat. no. Art. Nr.		PU VPE
	1-3718		1

driver for nail
Einschläger für Nagel



	Cat. no. Art. Nr.		PU VPE
	1-37011		1

slide hammer for humeral nail, cannulated
Gleithammer für Humerusnagel, kanüliert




	Cat. no. Art. Nr.		PU VPE
	1-1890	optional	1

nail length gauge
Messlehre



	Cat. no. Art. Nr.		PU VPE
	1-1098		1

diamond point awl
Ahle





TREU-Instrumente GmbH
take-off Gewerbepark 130 - 132
78579 Neuhausen o. E./Germany
Phone +49 (0) 7467 / 947 673 - 0
Fax +49 (0) 7467 / 947 673 - 50
info@treu.com
www.treu.com