





YOUR RELIABLE PARTNER IN TRAUMA CARE

WE FOCUS THE FUTURE

## LARGE FRAGMENT SYSTEM







Set Cat. No. T 340


Cat. no. Art. Nr.	included items	QTY St.
T 341	Tray for Basic Instrument Set	1
9-020	Drill Bit AO Shaft Ø 3.2mm x 145/120mm	3
9-026	Drill Bit AO Shaft Ø 4.5mm x 145/120mm	2
9-104	Countersink for Malleolar Screws Ø 3.2mm	1
9-092	Tap Handle for AO Shaft	1
9-078	Tap AO Shaft long thread Ø 4.5mm	2
9-079	Tap AO Shaft Ø 6.5mm	1
9-256	Tap Sleeve Ø 3.5mm	1
9-280	Insert Drill Sleeve 4.5mm / 3.2mm	1
9-260	Tap and Drill Sleeve Ø 4.5mm	1
9-144	Hexagonal Screwdriver Ø 3.5mm	1
9-176	Hexagonal Screwdriver AO Shaft Ø 3.5mm	1
9-114	Depth Gauge for Large Screws	1
9-116	Sharp Hook	1
9-262	Drill Guide for Round-Hole-Plate Ø 3.2mm	1
9-440	Tension Device 20mm	1
9-432	Socket Wrench width 11mm	1
9-1861	Combination Wrench 11mm	1
9-272	Neutral and Load Drill Guide Ø 3.2mm	1
3-070	Bending Templates 5 Holes standard 9cm	1
3-071	Bending Templates 7 Holes standard 12cm	1
3-072	Bending Templates 9 holes standard 15cm	1
9-106	Countersink for Large Screws	1
7-264	K-wire round 2.0mm dia x 150mm	10

Basic Instrument Set for Screws and Plates  
Basis Instrumentensatz für Schrauben und Platten



Sterilization Container (recommended)

15-1870	Bottom non-perforated 580 x 280 x 100 mm
15-1850	Lid perforated 580 x 280mm silver






## Set Cat. No. T 300

Cat. no. Art. Nr.	included items	QTY St.
T 111	Aluminium Case	1
T 301	Upper Tray	1
T 402	Lower Tray	1
9-020	Drill Bit AO Shaft Ø 3.2mm x 145/120mm	3
9-026	Drill Bit AO Shaft Ø 4.5mm x 145/120mm	2
9-104	Countersink for Malleolar Screws Ø 3.2mm	1
9-092	Tap Handle for AO Shaft	1
9-077	Tap AO Shaft short thread Ø 4.5mm	1
9-078	Tap AO Shaft long thread Ø 4.5mm	1
9-079	Tap AO Shaft Ø 6.5mm	1
9-256	Tap Sleeve Ø 3.5mm	1
9-280	Insert Drill Sleeve 4.5mm / 3.2mm	1
9-260	Tap and Drill Sleeve Ø 4.5mm	1
9-258	Pointed Drill Guide	1
9-144	Hexagonal Screwdriver Ø 3.5mm	1
9-176	Hexagonal Screwdriver AO Shaft Ø 3.5mm	1
9-114	Depth Gauge for Large Screws	1
9-116	Sharp Hook	1
9-430	Drill Sleeve for Tension Device	1
9-262	Drill Guide for Round-Hole-Plate Ø 3.2mm	1
9-436	Tension Device 8mm	1
9-440	Tension Device 20mm	1
9-432	Socket Wrench width 11mm	1
9-1861	Combination Wrench 11mm	1
9-272	Neutral and Load Drill Guide Ø 3.2mm	1
3-070	Bending Templates 5 Holes standard 9cm	1
3-071	Bending Templates 7 Holes standard 12cm	1
3-072	Bending Templates 9 Holes standard 15cm	1
9-106	Countersink for Large Screws	1
15-3351	Wrapping Drape cotton green 72 x 43cm	1

Basic Instrument Set for Screws and Plates  
Basis Instrumentensatz für Schrauben und Platten





Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-020	3.2mm	145mm	1

Drill Bit AO Shaft Ø 3.2mm x 145/120mm  
Bohrer AO Schaft Ø 3.2mm x 145/120mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-026	4.5mm	145mm	1

Drill Bit AO Shaft Ø 4.5mm x 145/120mm  
Bohrer AO Schaft Ø 4.5mm x 145/120mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-079	6.5mm	200mm	1

Tap AO Shaft Ø 6.5mm  
Gewindeschneider AO Schaft Ø 6.5mm




Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-078	4.5mm	130mm	1

Tap AO Shaft long thread Ø 4.5mm  
Gewindeschneider AO Schaft langes Gewinde Ø 4.5mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-256	3.5mm	-	1

Tap Sleeve Ø 3.5mm  
Gewebeschtzhülse Ø 3.5mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-260	4.5mm	-	1

Tap and Drill Sleeve Ø 4.5mm  
Gewebeschutzhülse und Bohrbuchse Ø 4.5mm




Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-262	3.2mm	-	1

Drill Guide for Round-Hole-Plate Ø 3.2mm  
Gewebeschutzhülse für Rundlochplatten Ø 3.2mm




Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-092	-	80mm	1

Tap Handle AO Shaft  
Vorbohrergriff AO Schaft



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-106	4.5mm	180mm	1

Countersink for Large Screws  
Kopfraumfräser für große Schrauben





Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-272	3.2mm	170mm	1

Neutral and Load Drill Guide Ø 3.2mm  
Neutral- und Kompressionsbohrbüchse Ø 3.2mm




Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-280	3.2mm	60mm	1

Insert Drill Sleeve 4.5mm / 3.2mm  
Steck-Bohrbüchse 4.5mm / 3.2mm




Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-258	-	-	1

Pointed Drill Guide  
Bohrbüchse mit Spitze



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-432	11mm	195mm	1

Socket Wrench 11mm  
Kardanschlüssel 11mm





Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-104	3.2mm	115mm	1

Countersink for Malleolar Screws Ø 3.2mm  
Kopfraumfräser für Malleolarschrauben Ø 3.2mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-144	3.5mm	200mm	1

Hexagonal Screwdriver Ø 3.5mm  
Hexagonaler Schraubendreher Ø 3.5mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-176	3.5mm	100mm	1

Hexagonal Screwdriver AO Shaft Ø 3.5mm  
Schraubenzieherschaft AO Ø 3.5mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-430	-	-	1

Drill Sleeve for Tension Device  
Bohrlehre für Plattenspanner



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-1861	11mm	160mm	1

Combination Wrench 11mm  
Gabelringschlüssel 11mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-114	-	-	1

Depth Gauge for Large Screws  
Schraubenmesser für große Schrauben



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-116	-	150mm	1

Sharp Hook  
Scharfer Haken




Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
4-184	-	-	1

Screw Holding Forceps  
Schraubenhaltezange




Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-070	5	90mm	1
3-071	7	120mm	1
3-072	9	150mm	1

Bending Templates standard  
Biegeschablone Standard



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
9-436	8mm	60mm	1

Tension Device 8mm  
Plattenspanner 8mm




60mm




			
Cat. no. Art. Nr.	∅	Length Länge	PU VPE
9-440	20mm	-	1


Tension Device 20mm  
Plattenspanner 20mm



			
Cat. no. Art. Nr.	∅	Length Länge	PU VPE
9-416	-	-	1

Bending Iron for Standard Plates  
Biegeisen für Standardplatten



			
Cat. no. Art. Nr.	∅	Length Länge	PU VPE
9-442	-	-	1

Pin Wrench for Tension Device and Distractor  
Schlüssel für Plattenspanner und Distraktor




			
Cat. no. Art. Nr.	∅	Length Länge	PU VPE
9-420	-	-	1

Plate Bending Press  
Biegepresse





Set Cat. No. T 336

Screw Set in Graphic Case  
Schraubensatz im Graphic Case



Sterilization Container (recommended)	
15-1952	Bottom non-perforated 285 x 280 x 135 mm
15-1930	Lid perforated 285 x 280mm silver



 Set Cat. No. T 336		
Cat. no. Art. Nr.	included items	QTY St.
T 337	Tray for Screw Set	1
<b>4.5mm Cortical Screws</b>		
3-200	Cortical Screw 4.5 x 14mm	6
3-202	Cortical Screw 4.5 x 16mm	6
3-204	Cortical Screw 4.5 x 18mm	6
3-206	Cortical Screw 4.5 x 20mm	6
3-208	Cortical Screw 4.5 x 22mm	6
3-210	Cortical Screw 4.5 x 24mm	6
3-212	Cortical Screw 4.5 x 26mm	7
3-214	Cortical Screw 4.5 x 28mm	7
3-216	Cortical Screw 4.5 x 30mm	14
3-218	Cortical Screw 4.5 x 32mm	14
3-220	Cortical Screw 4.5 x 34mm	14
3-222	Cortical Screw 4.5 x 36mm	14
3-224	Cortical Screw 4.5 x 38mm	7
3-226	Cortical Screw 4.5 x 40mm	14
3-228	Cortical Screw 4.5 x 42mm	7
3-230	Cortical Screw 4.5 x 44mm	7
3-232	Cortical Screw 4.5 x 46mm	7
3-234	Cortical Screw 4.5 x 48mm	3
3-238	Cortical Screw 4.5 x 52mm	3
3-242	Cortical Screw 4.5 x 56mm	3
3-246	Cortical Screw 4.5 x 60mm	3
3-248	Cortical Screw 4.5 x 65mm	3
3-250	Cortical Screw 4.5 x 70mm	3
<b>6.5mm / 16mm Cancellous Screws</b>		
3-382	Cancellous Screw 6.5 x 16mm thr. 30mm	2
3-384	Cancellous Screw 6.5 x 16mm thr. 35mm	2
3-386	Cancellous Screw 6.5 x 16mm thr. 40mm	2
3-388	Cancellous Screw 6.5 x 16mm thr. 45mm	2
3-390	Cancellous Screw 6.5 x 16mm thr. 50mm	2
3-392	Cancellous Screw 6.5 x 16mm thr. 55mm	2
3-394	Cancellous Screw 6.5 x 16mm thr. 60mm	2
3-396	Cancellous Screw 6.5 x 16mm thr. 65mm	2
3-398	Cancellous Screw 6.5 x 16mm thr. 70mm	2
3-400	Cancellous Screw 6.5 x 75mm thr. 16mm	2
3-402	Cancellous Screw 6.5 x 80mm thr. 16mm	2
3-404	Cancellous Screw 6.5 x 85mm thr. 16mm	2
3-406	Cancellous Screw 6.5 x 90mm thr. 16mm	2
3-408	Cancellous Screw 6.5 x 95mm thr. 16mm	2
3-410	Cancellous Screw 6.5 x 100mm thr. 16mm	2
3-412	Cancellous Screw 6.5 x 105mm thr. 16mm	2
3-414	Cancellous Screw 6.5 x 110mm thr. 16mm	2

Cat. no. Art. Nr.	included items	QTY St.
<b>6.5mm / 32mm Cancellous Screws</b>		
3-436	Cancellous Screw 6.5 x 45mm thr. 32mm	3
3-438	Cancellous Screw 6.5 x 50mm thr. 32mm	3
3-440	Cancellous Screw 6.5 x 55mm thr. 32mm	3
3-442	Cancellous Screw 6.5 x 60mm thr. 32mm	3
3-444	Cancellous Screw 6.5 x 65mm thr. 32mm	3
3-446	Cancellous Screw 6.5 x 70mm thr. 32mm	3
3-448	Cancellous Screw 6.5 x 75mm thr. 32mm	3
3-450	Cancellous Screw 6.5 x 80mm thr. 32mm	3
3-452	Cancellous Screw 6.5 x 85mm thr. 32mm	3
3-454	Cancellous Screw 6.5 x 90mm thr. 32mm	3
3-456	Cancellous Screw 6.5 x 95mm thr. 32mm	3
3-458	Cancellous Screw 6.5 x 32mm thr. 100mm	3
3-460	Cancellous Screw 6.5 x 32mm thr. 105mm	3
3-462	Cancellous Screw 6.5 x 32mm thr. 110mm	3
<b>4.5mm Malleolar Screws</b>		
3-352	Malleolar Screw 4.5 x 25mm	3
3-354	Malleolar Screw 4.5 x 30mm	3
3-356	Malleolar Screw 4.5 x 35mm	3
3-358	Malleolar Screw 4.5 x 40mm	3
3-360	Malleolar Screw 4.5 x 45mm	3
3-362	Malleolar Screw 4.5 x 50mm	3
3-364	Malleolar Screw 4.5 x 55mm	3
3-366	Malleolar Screw 4.5 x 60mm	3
3-368	Malleolar Screw 4.5 x 65mm	3
3-370	Malleolar Screw 4.5 x 70mm	3
<b>6.5mm Cancellous Screws full thread</b>		
3-480	Cancellous Screw 6.5 x 25mm full thr.	4
3-482	Cancellous Screw 6.5 x 30mm full thr.	4
3-484	Cancellous Screw 6.5 x 35mm full thr.	4
3-486	Cancellous Screw 6.5 x 40mm full thr.	4
3-488	Cancellous Screw 6.5 x 45mm full thr.	2
3-490	Cancellous Screw 6.5 x 50mm full thr.	2
3-492	Cancellous Screw 6.5 x 55mm full thr.	2
3-494	Cancellous Screw 6.5 x 60mm full thr.	2
3-496	Cancellous Screw 6.5 x 65mm full thr.	2
3-498	Cancellous Screw 6.5 x 70mm full thr.	2
3-500	Cancellous Screw 6.5 x 75mm full thr.	2
3-502	Cancellous Screw 6.5 x 80mm full thr.	2
3-504	Cancellous Screw 6.5 x 85mm full thr.	2
3-506	Cancellous Screw 6.5 x 90mm full thr.	2
3-508	Cancellous Screw 6.5 x 95mm full thr.	2
3-512	Cancellous Screw 6.5 x 105mm full thr.	2
3-514	Cancellous Screw 6.5 x 110mm full thr.	2
<b>Equipment</b>		
3-522	Washer for Cancellous ext. Ø 13mm	6
4-184	Screw Holding Forceps	1




Set Cat. No. T 3360

Titanium Screw Set in Graphic Case  
Titan Schraubensatz im Graphic Case



Sterilization Container (recommended)	
15-1952	Bottom non-perforated 285 x 280 x 135 mm
15-1930	Lid perforated 285 x 280mm silver



 Set Cat. No. T 3360		
Cat. no. Art. Nr.	included items	QTY St.
T 337	Tray for Screw Set	1
4.5mm Cortical Screws		
3-20000	Cortical Screw 4.5 x 14mm Titan	6
3-20200	Cortical Screw 4.5 x 16mm Titan	6
3-20400	Cortical Screw 4.5 x 18mm Titan	6
3-20600	Cortical Screw 4.5 x 20mm Titan	6
3-20800	Cortical Screw 4.5 x 22mm Titan	6
3-21000	Cortical Screw 4.5 x 24mm Titan	6
3-21200	Cortical Screw 4.5 x 26mm Titan	7
3-21400	Cortical Screw 4.5 x 28mm Titan	7
3-21600	Cortical Screw 4.5 x 30mm Titan	14
3-21800	Cortical Screw 4.5 x 32mm Titan	14
3-22000	Cortical Screw 4.5 x 34mm Titan	14
3-22200	Cortical Screw 4.5 x 36mm Titan	14
3-22400	Cortical Screw 4.5 x 38mm Titan	7
3-22600	Cortical Screw 4.5 x 40mm Titan	14
3-22800	Cortical Screw 4.5 x 42mm Titan	7
3-23000	Cortical Screw 4.5 x 44mm Titan	7
3-23200	Cortical Screw 4.5 x 46mm Titan	7
3-23400	Cortical Screw 4.5 x 48mm Titan	3
3-23800	Cortical Screw 4.5 x 52mm Titan	3
3-24200	Cortical Screw 4.5 x 56mm Titan	3
3-24600	Cortical Screw 4.5 x 60mm Titan	3
3-24800	Cortical Screw 4.5 x 65mm Titan	3
3-25000	Cortical Screw 4.5 x 70mm Titan	3
6.5mm / 16mm Cancellous Screws		
3-38200	Cancellous Screw 6.5 x 16mm thr. 30mm Titan	2
3-38400	Cancellous Screw 6.5 x 16mm thr. 35mm Titan	2
3-38600	Cancellous Screw 6.5 x 16mm thr. 40mm Titan	2
3-38800	Cancellous Screw 6.5 x 16mm thr. 45mm Titan	2
3-39000	Cancellous Screw 6.5 x 16mm thr. 50mm Titan	2
3-39200	Cancellous Screw 6.5 x 16mm thr. 55mm Titan	2
3-39400	Cancellous Screw 6.5 x 16mm thr. 60mm Titan	2
3-39600	Cancellous Screw 6.5 x 16mm thr. 65mm Titan	2
3-39800	Cancellous Screw 6.5 x 16mm thr. 70mm Titan	2
3-40000	Cancellous Screw 6.5 x 75mm thr. 16mm Titan	2
3-40200	Cancellous Screw 6.5 x 80mm thr. 16mm Titan	2
3-40400	Cancellous Screw 6.5 x 85mm thr. 16mm Titan	2
3-40600	Cancellous Screw 6.5 x 90mm thr. 16mm Titan	2
3-40800	Cancellous Screw 6.5 x 95mm thr. 16mm Titan	2
3-41000	Cancellous Screw 6.5 x 100mm thr. 16mm Titan	2
3-41200	Cancellous Screw 6.5 x 105mm thr. 16mm Titan	2
3-41400	Cancellous Screw 6.5 x 110mm thr. 16mm Titan	2

Cat. no. Art. Nr.	included items	QTY St.
6.5mm / 32mm Cancellous Screws		
3-43600	Cancellous Screw 6.5 x 45mm thr. 32mm Titan	3
3-43800	Cancellous Screw 6.5 x 50mm thr. 32mm Titan	3
3-44000	Cancellous Screw 6.5 x 55mm thr. 32mm Titan	3
3-44200	Cancellous Screw 6.5 x 60mm thr. 32mm Titan	3
3-44400	Cancellous Screw 6.5 x 65mm thr. 32mm Titan	3
3-44600	Cancellous Screw 6.5 x 70mm thr. 32mm Titan	3
3-44800	Cancellous Screw 6.5 x 75mm thr. 32mm Titan	3
3-45000	Cancellous Screw 6.5 x 80mm thr. 32mm Titan	3
3-45200	Cancellous Screw 6.5 x 85mm thr. 32mm Titan	3
3-45400	Cancellous Screw 6.5 x 90mm thr. 32mm Titan	3
3-45600	Cancellous Screw 6.5 x 95mm thr. 32mm Titan	3
3-45800	Cancellous Screw 6.5 x 32mm thr. 100mm Titan	3
3-46000	Cancellous Screw 6.5 x 32mm thr. 105mm Titan	3
3-46200	Cancellous Screw 6.5 x 32mm thr. 110mm Titan	3
4.5mm Shank Screws		
3-28600	Shank Screw 4.5 x 28mm Titan	2
3-28800	Shank Screw 4.5 x 30mm Titan	2
3-29000	Shank Screw 4.5 x 32mm Titan	2
3-29200	Shank Screw 4.5 x 34mm Titan	2
3-29400	Shank Screw 4.5 x 36mm Titan	2
3-29600	Shank Screw 4.5 x 38mm Titan	2
3-29800	Shank Screw 4.5 x 40mm Titan	2
3-30000	Shank Screw 4.5 x 42mm Titan	2
3-30200	Shank Screw 4.5 x 44mm Titan	2
3-30400	Shank Screw 4.5 x 46mm Titan	2
6.5mm Cancellous Screws full thread		
3-48000	Cancellous Screw 6.5 x 25mm full thr. Titan	4
3-48200	Cancellous Screw 6.5 x 30mm full thr. Titan	4
3-48400	Cancellous Screw 6.5 x 35mm full thr. Titan	4
3-48600	Cancellous Screw 6.5 x 40mm full thr. Titan	4
3-48800	Cancellous Screw 6.5 x 45mm full thr. Titan	2
3-49000	Cancellous Screw 6.5 x 50mm full thr. Titan	2
3-49200	Cancellous Screw 6.5 x 55mm full thr. Titan	2
3-49400	Cancellous Screw 6.5 x 60mm full thr. Titan	2
3-49600	Cancellous Screw 6.5 x 65mm full thr. Titan	2
3-49800	Cancellous Screw 6.5 x 70mm full thr. Titan	2
3-50000	Cancellous Screw 6.5 x 75mm full thr. Titan	2
3-50200	Cancellous Screw 6.5 x 80mm full thr. Titan	2
3-50400	Cancellous Screw 6.5 x 85mm full thr. Titan	2
3-50600	Cancellous Screw 6.5 x 90mm full thr. Titan	2
3-50800	Cancellous Screw 6.5 x 95mm full thr. Titan	2
3-51200	Cancellous Screw 6.5 x 105mm full thr. Titan	2
3-51400	Cancellous Screw 6.5 x 110mm full thr. Titan	2
Equipment		
3-5220	Washer Titanium for Cancellous ext. Ø 13mm	6
4-184	Screw Holding Forceps	1




Set Cat. No. T 306

Screw Set  
Schraubenset





 Set Cat. No. T 306		
Cat. no. Art. Nr.	included items	QTY St.
T 111	Aluminium Case red	1
T 307	Tray (Aluminium)	1
T 308	Tray (Plastic)	1
<b>4.5mm Cortical Screws</b>		
3-200	Cortical Screw 4.5 x 14mm	7
3-202	Cortical Screw 4.5 x 16mm	8
3-204	Cortical Screw 4.5 x 18mm	8
3-206	Cortical Screw 4.5 x 20mm	8
3-208	Cortical Screw 4.5 x 22mm	8
3-210	Cortical Screw 4.5 x 24mm	8
3-212	Cortical Screw 4.5 x 26mm	8
3-214	Cortical Screw 4.5 x 28mm	8
3-216	Cortical Screw 4.5 x 30mm	16
3-218	Cortical Screw 4.5 x 32mm	12
3-220	Cortical Screw 4.5 x 34mm	12
3-222	Cortical Screw 4.5 x 36mm	12
3-224	Cortical Screw 4.5 x 38mm	12
3-226	Cortical Screw 4.5 x 40mm	12
3-228	Cortical Screw 4.5 x 42mm	6
3-230	Cortical Screw 4.5 x 44mm	12
3-232	Cortical Screw 4.5 x 46mm	6
3-234	Cortical Screw 4.5 x 48mm	6
3-238	Cortical Screw 4.5 x 52mm	6
3-242	Cortical Screw 4.5 x 56mm	6
3-246	Cortical Screw 4.5 x 60mm	5
3-248	Cortical Screw 4.5 x 65mm	4
3-250	Cortical Screw 4.5 x 70mm	4
<b>6.5mm / 16mm Cancellous Screws</b>		
3-382	Cancellous Screw 6.5 x 16mm thr. 30mm	2
3-384	Cancellous Screw 6.5 x 16mm thr. 35mm	2
3-386	Cancellous Screw 6.5 x 16mm thr. 40mm	2
3-388	Cancellous Screw 6.5 x 16mm thr. 45mm	2
3-390	Cancellous Screw 6.5 x 16mm thr. 50mm	2
3-392	Cancellous Screw 6.5 x 16mm thr. 55mm	2
3-394	Cancellous Screw 6.5 x 16mm thr. 60mm	2
3-396	Cancellous Screw 6.5 x 16mm thr. 65mm	3
3-398	Cancellous Screw 6.5 x 16mm thr. 70mm	3
3-400	Cancellous Screw 6.5 x 75mm thr. 16mm	3
3-402	Cancellous Screw 6.5 x 80mm thr. 16mm	3
3-404	Cancellous Screw 6.5 x 85mm thr. 16mm	3
3-406	Cancellous Screw 6.5 x 90mm thr. 16mm	3
3-408	Cancellous Screw 6.5 x 95mm thr. 16mm	3
3-410	Cancellous Screw 6.5 x 100mm thr. 16mm	3
3-412	Cancellous Screw 6.5 x 105mm thr. 16mm	3
3-414	Cancellous Screw 6.5 x 110mm thr. 16mm	3

Cat. no. Art. Nr.	included items	QTY St.
<b>6.5mm / 32mm Cancellous Screws</b>		
3-436	Cancellous Screw 6.5 x 45mm thr. 32mm	4
3-438	Cancellous Screw 6.5 x 50mm thr. 32mm	4
3-440	Cancellous Screw 6.5 x 55mm thr. 32mm	4
3-442	Cancellous Screw 6.5 x 60mm thr. 32mm	4
3-444	Cancellous Screw 6.5 x 65mm thr. 32mm	3
3-446	Cancellous Screw 6.5 x 70mm thr. 32mm	3
3-450	Cancellous Screw 6.5 x 80mm thr. 32mm	3
3-452	Cancellous Screw 6.5 x 85mm thr. 32mm	3
3-454	Cancellous Screw 6.5 x 90mm thr. 32mm	3
<b>6.5mm Cancellous Screws full thread</b>		
3-480	Cancellous Screw 6.5 x 25mm full thr.	6
3-482	Cancellous Screw 6.5 x 30mm full thr.	4
3-484	Cancellous Screw 6.5 x 35mm full thr.	4
3-486	Cancellous Screw 6.5 x 40mm full thr.	4
<b>Equipment</b>		
4-184	Screw Holding Forceps	1
3-522	Washer for Cancellous ext. Ø 13mm	6
15-3351	Wrapping Drape cotton green 72 x 43cm	1



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-200	4.5mm	14mm	1
3-202	4.5mm	16mm	1
3-204	4.5mm	18mm	1
3-206	4.5mm	20mm	1
3-208	4.5mm	22mm	1
3-210	4.5mm	24mm	1
3-212	4.5mm	26mm	1
3-214	4.5mm	28mm	1
3-216	4.5mm	30mm	1
3-218	4.5mm	32mm	1
3-220	4.5mm	34mm	1
3-222	4.5mm	36mm	1
3-224	4.5mm	38mm	1
3-226	4.5mm	40mm	1
3-228	4.5mm	42mm	1
3-230	4.5mm	44mm	1
3-232	4.5mm	46mm	1
3-234	4.5mm	48mm	1
3-236	4.5mm	50mm	1
3-238	4.5mm	52mm	1
3-240	4.5mm	54mm	1
3-242	4.5mm	56mm	1
3-244	4.5mm	58mm	1
3-246	4.5mm	60mm	1
3-248	4.5mm	65mm	1
3-250	4.5mm	70mm	1
3-254	4.5mm	75mm	1
3-260	4.5mm	80mm	1
3-262	4.5mm	85mm	1
3-264	4.5mm	90mm	1
3-266	4.5mm	95mm	1
3-268	4.5mm	100mm	1
3-269	4.5mm	105mm	1
3-270	4.5mm	110mm	1
3-271	4.5mm	115mm	1
3-272	4.5mm	120mm	1

### Cortical Screw 4.5mm (Stainless Steel) Kortikalisschraube 4.5mm (Stahl)

Screw Diameter • Schraubendurchmesser	4.50mm
Core Diameter • Kerndurchmesser	3.00mm
Head Diameter • Kopfdurchmesser	8.00mm
Drill Bit for threaded Hole • Bohrer für Kernloch	3.20mm





Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-20000	4.5mm	14mm	1
3-20200	4.5mm	16mm	1
3-20400	4.5mm	18mm	1
3-20600	4.5mm	20mm	1
3-20800	4.5mm	22mm	1
3-21000	4.5mm	24mm	1
3-21200	4.5mm	26mm	1
3-21400	4.5mm	28mm	1
3-21600	4.5mm	30mm	1
3-21800	4.5mm	32mm	1
3-22000	4.5mm	34mm	1
3-22200	4.5mm	36mm	1
3-22400	4.5mm	38mm	1
3-22600	4.5mm	40mm	1
3-22800	4.5mm	42mm	1
3-23000	4.5mm	44mm	1
3-23200	4.5mm	46mm	1
3-23400	4.5mm	48mm	1
3-23600	4.5mm	50mm	1
3-23800	4.5mm	52mm	1
3-24000	4.5mm	54mm	1
3-24200	4.5mm	56mm	1
3-24400	4.5mm	58mm	1
3-24600	4.5mm	60mm	1
3-24800	4.5mm	65mm	1
3-25000	4.5mm	70mm	1
3-25400	4.5mm	75mm	1
3-26000	4.5mm	80mm	1
3-26200	4.5mm	85mm	1
3-26400	4.5mm	90mm	1
3-26600	4.5mm	95mm	1
3-26800	4.5mm	100mm	1
3-26900	4.5mm	105mm	1
3-27000	4.5mm	110mm	1
3-27100	4.5mm	115mm	1
3-27200	4.5mm	120mm	1

### Cortical Screw 4.5mm (Titanium) Kortikalisschraube 4.5mm (Titan)

Screw Diameter • Schraubendurchmesser	4.50mm
Core Diameter • Kerndurchmesser	3.00mm
Head Diameter • Kopfdurchmesser	8.00mm
Drill Bit for threaded Hole • Bohrer für Kernloch	3.20mm






Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-1200	4.5mm	14mm	1
3-1202	4.5mm	16mm	1
3-1204	4.5mm	18mm	1
3-1206	4.5mm	20mm	1
3-1208	4.5mm	22mm	1
3-1210	4.5mm	24mm	1
3-1212	4.5mm	26mm	1
3-1214	4.5mm	28mm	1
3-1216	4.5mm	30mm	1
3-1218	4.5mm	32mm	1
3-1220	4.5mm	34mm	1
3-1222	4.5mm	36mm	1
3-1224	4.5mm	38mm	1
3-1226	4.5mm	40mm	1
3-1228	4.5mm	42mm	1
3-1230	4.5mm	44mm	1
3-1232	4.5mm	46mm	1
3-1234	4.5mm	48mm	1
3-1236	4.5mm	50mm	1
3-1238	4.5mm	52mm	1
3-1240	4.5mm	54mm	1
3-1242	4.5mm	56mm	1
3-1244	4.5mm	58mm	1
3-1246	4.5mm	60mm	1
3-1248	4.5mm	65mm	1
3-1250	4.5mm	70mm	1
3-1252	4.5mm	75mm	1
3-1254	4.5mm	80mm	1
3-1256	4.5mm	85mm	1
3-1258	4.5mm	90mm	1
3-1260	4.5mm	95mm	1
3-1262	4.5mm	100mm	1
3-1264	4.5mm	105mm	1
3-1266	4.5mm	110mm	1
3-1268	4.5mm	115mm	1
3-1270	4.5mm	120mm	1

### Cortical Screw 4.5mm self-cutting (Stainless Steel) Kortikalisschraube 4.5mm selbstschneidend (Stahl)

Screw Diameter • Schraubendurchmesser	4.50mm
Core Diameter • Kerndurchmesser	3.00mm
Head Diameter • Kopfdurchmesser	8.00mm
Drill Bit for threaded Hole • Bohrer für Kernloch	3.20mm




self-cutting Screws, tapping not needed



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-120000	4.5mm	14mm	1
3-120200	4.5mm	16mm	1
3-120400	4.5mm	18mm	1
3-120600	4.5mm	20mm	1
3-120800	4.5mm	22mm	1
3-121000	4.5mm	24mm	1
3-121200	4.5mm	26mm	1
3-121400	4.5mm	28mm	1
3-121600	4.5mm	30mm	1
3-121800	4.5mm	32mm	1
3-122000	4.5mm	34mm	1
3-122200	4.5mm	36mm	1
3-122400	4.5mm	38mm	1
3-122600	4.5mm	40mm	1
3-122800	4.5mm	42mm	1
3-123000	4.5mm	44mm	1
3-123200	4.5mm	46mm	1
3-123400	4.5mm	48mm	1
3-123600	4.5mm	50mm	1
3-123800	4.5mm	52mm	1
3-124000	4.5mm	54mm	1
3-124200	4.5mm	56mm	1
3-124400	4.5mm	58mm	1
3-124600	4.5mm	60mm	1
3-124800	4.5mm	65mm	1
3-125000	4.5mm	70mm	1
3-125200	4.5mm	75mm	1
3-125400	4.5mm	80mm	1
3-125600	4.5mm	85mm	1
3-125800	4.5mm	90mm	1
3-126000	4.5mm	95mm	1
3-126200	4.5mm	100mm	1
3-126400	4.5mm	105mm	1
3-126600	4.5mm	110mm	1
3-126800	4.5mm	115mm	1
3-127000	4.5mm	120mm	1

### Cortical Screw 4.5mm self-cutting (Titanium) Kortikalisschraube 4.5mm selbstschneidend (Titan)

Screw Diameter • Schraubendurchmesser	4.50mm
Core Diameter • Kerndurchmesser	3.00mm
Head Diameter • Kopfdurchmesser	8.00mm
Drill Bit for threaded Hole • Bohrer für Kernloch	3.20mm



self-cutting Screws, tapping not needed



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-280	4.5mm	22mm	1
3-282	4.5mm	24mm	1
3-284	4.5mm	26mm	1
3-286	4.5mm	28mm	1
3-288	4.5mm	30mm	1
3-290	4.5mm	32mm	1
3-292	4.5mm	34mm	1
3-294	4.5mm	36mm	1
3-296	4.5mm	38mm	1
3-298	4.5mm	40mm	1
3-300	4.5mm	42mm	1
3-302	4.5mm	44mm	1
3-304	4.5mm	46mm	1
3-306	4.5mm	48mm	1
3-308	4.5mm	50mm	1
3-310	4.5mm	52mm	1
3-312	4.5mm	54mm	1
3-314	4.5mm	56mm	1
3-316	4.5mm	58mm	1
3-318	4.5mm	60mm	1
3-320	4.5mm	65mm	1
3-322	4.5mm	70mm	1
3-324	4.5mm	75mm	1
3-326	4.5mm	80mm	1
3-328	4.5mm	85mm	1
3-330	4.5mm	90mm	1
3-332	4.5mm	95mm	1
3-334	4.5mm	100mm	1
3-336	4.5mm	105mm	1
3-338	4.5mm	110mm	1

## Shank Screw 4.5mm (Stainless Steel) Schafschraube 4.5mm (Stahl)

Screw Diameter • Schraubendurchmesser	4.50mm
Core Diameter • Kerndurchmesser	3.10mm
Shaft Diameter • Schaftdurchmesser	4.50mm
Head Diameter • Kopfdurchmesser	8.00mm
Drill Bit for threaded Hole • Bohrer für Kernloch	3.20mm





Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-28000	4.5mm	22mm	1
3-28200	4.5mm	24mm	1
3-28400	4.5mm	26mm	1
3-28600	4.5mm	28mm	1
3-28800	4.5mm	30mm	1
3-29000	4.5mm	32mm	1
3-29200	4.5mm	34mm	1
3-29400	4.5mm	36mm	1
3-29600	4.5mm	38mm	1
3-29800	4.5mm	40mm	1
3-30000	4.5mm	42mm	1
3-30200	4.5mm	44mm	1
3-30400	4.5mm	46mm	1
3-30600	4.5mm	48mm	1
3-30800	4.5mm	50mm	1
3-31000	4.5mm	52mm	1
3-31200	4.5mm	54mm	1
3-31400	4.5mm	56mm	1
3-31600	4.5mm	58mm	1
3-31800	4.5mm	60mm	1
3-32000	4.5mm	65mm	1
3-32200	4.5mm	70mm	1
3-32400	4.5mm	75mm	1
3-32600	4.5mm	80mm	1
3-32800	4.5mm	85mm	1
3-33000	4.5mm	90mm	1
3-33200	4.5mm	95mm	1
3-33400	4.5mm	100mm	1
3-33600	4.5mm	105mm	1
3-33800	4.5mm	110mm	1

### Shank Screw 4.5mm (Titanium) Schaftschraube 4.5mm (Titan)


Screw Diameter • Schraubendurchmesser	4.50mm
Core Diameter • Kerndurchmesser	3.10mm
Shaft Diameter • Schaftdurchmesser	4.50mm
Head Diameter • Kopfdurchmesser	8.00mm
Drill Bit for threaded Hole • Bohrer für Kernloch	3.20mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-1280	4.5mm	22mm	1
3-1282	4.5mm	24mm	1
3-1284	4.5mm	26mm	1
3-1286	4.5mm	28mm	1
3-1288	4.5mm	30mm	1
3-1290	4.5mm	32mm	1
3-1292	4.5mm	34mm	1
3-1294	4.5mm	36mm	1
3-1296	4.5mm	38mm	1
3-1298	4.5mm	40mm	1
3-1300	4.5mm	42mm	1
3-1302	4.5mm	44mm	1
3-1304	4.5mm	46mm	1
3-1306	4.5mm	48mm	1
3-1308	4.5mm	50mm	1
3-1310	4.5mm	52mm	1
3-1312	4.5mm	54mm	1
3-1314	4.5mm	56mm	1
3-1316	4.5mm	58mm	1
3-1318	4.5mm	60mm	1
3-1320	4.5mm	65mm	1
3-1322	4.5mm	70mm	1
3-1324	4.5mm	75mm	1
3-1326	4.5mm	80mm	1
3-1328	4.5mm	85mm	1
3-1330	4.5mm	90mm	1
3-1332	4.5mm	95mm	1
3-1334	4.5mm	100mm	1
3-1336	4.5mm	105mm	1
3-1338	4.5mm	110mm	1

Shank Screw 4.5mm self-cutting (Stainless Steel)  
 Schaftschraube 4.5mm selbstschneidend (Stahl)

Screw Diameter • Schraubendurchmesser 4.50mm  
 Core Diameter • Kerndurchmesser 3.10mm  
 Shaft Diameter • Schaftdurchmesser 4.50mm  
 Head Diameter • Kopfdurchmesser 8.00mm  
 Drill Bit for threaded Hole • Bohrer für Kernloch 3.20mm




self-cutting Screws, tapping not needed

Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-350	4.5mm	20mm	1
3-352	4.5mm	25mm	1
3-354	4.5mm	30mm	1
3-356	4.5mm	35mm	1
3-358	4.5mm	40mm	1
3-360	4.5mm	45mm	1
3-362	4.5mm	50mm	1
3-364	4.5mm	55mm	1
3-366	4.5mm	60mm	1
3-368	4.5mm	65mm	1
3-370	4.5mm	70mm	1
3-371	4.5mm	75mm	1
3-372	4.5mm	80mm	1

Malleolar Screw 4.5mm (Stainless Steel)  
 Malleolarschraube 4.5mm (Stahl)

Screw Diameter • Schraubendurchmesser 4.50mm  
 Core Diameter • Kerndurchmesser 3.00mm  
 Shaft Diameter • Schaftdurchmesser 3.00mm  
 Head Diameter • Kopfdurchmesser 8.00mm  
 Drill Bit for threaded Hole • Bohrer für Kernloch 3.20mm





Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-382	6.5mm	30mm	1
3-384	6.5mm	35mm	1
3-386	6.5mm	40mm	1
3-388	6.5mm	45mm	1
3-390	6.5mm	50mm	1
3-392	6.5mm	55mm	1
3-394	6.5mm	60mm	1
3-396	6.5mm	65mm	1
3-398	6.5mm	70mm	1
3-400	6.5mm	75mm	1
3-402	6.5mm	80mm	1
3-404	6.5mm	85mm	1
3-406	6.5mm	90mm	1
3-408	6.5mm	95mm	1
3-410	6.5mm	100mm	1
3-412	6.5mm	105mm	1
3-414	6.5mm	110mm	1
3-415	6.5mm	115mm	1
3-416	6.5mm	120mm	1

**Cancellous Screw 6.5mm (Stainless Steel)**  
Spongiosaschraube 6.5mm (Stahl)

Screw Diameter • Schraubendurchmesser 6.50mm  
 Core Diameter • Kerndurchmesser 3.00mm  
 Shaft Diameter • Schaftdurchmesser 4.50mm  
 Head Diameter • Kopfdurchmesser 8.00mm  
 Drill Bit for threaded Hole • Bohrer für Kernloch 3.20mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-38200	6.5mm	30mm	1
3-38400	6.5mm	35mm	1
3-38600	6.5mm	40mm	1
3-38800	6.5mm	45mm	1
3-39000	6.5mm	50mm	1
3-39200	6.5mm	55mm	1
3-39400	6.5mm	60mm	1
3-39600	6.5mm	65mm	1
3-39800	6.5mm	70mm	1
3-40000	6.5mm	75mm	1
3-40200	6.5mm	80mm	1
3-40400	6.5mm	85mm	1
3-40600	6.5mm	90mm	1
3-40800	6.5mm	95mm	1
3-41000	6.5mm	100mm	1
3-41200	6.5mm	105mm	1
3-41400	6.5mm	110mm	1
3-41500	6.5mm	115mm	1
3-41600	6.5mm	120mm	1

**Cancellous Screw 6.5mm (Titanium)**  
Spongiosaschraube 6.5mm (Titan)

Screw Diameter • Schraubendurchmesser 6.50mm  
 Core Diameter • Kerndurchmesser 3.00mm  
 Shaft Diameter • Schaftdurchmesser 4.50mm  
 Head Diameter • Kopfdurchmesser 8.00mm  
 Drill Bit for threaded Hole • Bohrer für Kernloch 3.20mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-434	6.5mm	40mm	1
3-436	6.5mm	45mm	1
3-438	6.5mm	50mm	1
3-440	6.5mm	55mm	1
3-442	6.5mm	60mm	1
3-444	6.5mm	65mm	1
3-446	6.5mm	70mm	1
3-448	6.5mm	75mm	1
3-450	6.5mm	80mm	1
3-452	6.5mm	85mm	1
3-454	6.5mm	90mm	1
3-456	6.5mm	95mm	1
3-458	6.5mm	100mm	1
3-460	6.5mm	105mm	1
3-462	6.5mm	110mm	1
3-464	6.5mm	115mm	1
3-466	6.5mm	120mm	1

**Cancellous Screw 6.5mm (Stainless Steel)**  
Spongiosaschraube 6.5mm (Stahl)

- Screw Diameter • Schraubendurchmesser 6.50mm
- Core Diameter • Kerndurchmesser 3.00mm
- Shaft Diameter • Schaftdurchmesser 4.50mm
- Head Diameter • Kopfdurchmesser 8.00mm
- Drill Bit for threaded Hole • Bohrer für Kernloch 3.20mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-43400	6.5mm	40mm	1
3-43600	6.5mm	45mm	1
3-43800	6.5mm	50mm	1
3-44000	6.5mm	55mm	1
3-44200	6.5mm	60mm	1
3-44400	6.5mm	65mm	1
3-44600	6.5mm	70mm	1
3-44800	6.5mm	75mm	1
3-45000	6.5mm	80mm	1
3-45200	6.5mm	85mm	1
3-45400	6.5mm	90mm	1
3-45600	6.5mm	95mm	1
3-45800	6.5mm	100mm	1
3-46000	6.5mm	105mm	1
3-46200	6.5mm	110mm	1
3-46400	6.5mm	115mm	1
3-46600	6.5mm	120mm	1

**Cancellous Screw 6.5mm (Titanium)**  
Spongiosaschraube 6.5mm (Titan)

- Screw Diameter • Schraubendurchmesser 6.50mm
- Core Diameter • Kerndurchmesser 3.00mm
- Shaft Diameter • Schaftdurchmesser 4.50mm
- Head Diameter • Kopfdurchmesser 8.00mm
- Drill Bit for threaded Hole • Bohrer für Kernloch 3.20mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-480	6.5mm	25mm	1
3-482	6.5mm	30mm	1
3-484	6.5mm	35mm	1
3-486	6.5mm	40mm	1
3-488	6.5mm	45mm	1
3-490	6.5mm	50mm	1
3-492	6.5mm	55mm	1
3-494	6.5mm	60mm	1
3-496	6.5mm	65mm	1
3-498	6.5mm	70mm	1
3-500	6.5mm	75mm	1
3-502	6.5mm	80mm	1
3-504	6.5mm	85mm	1
3-506	6.5mm	90mm	1
3-508	6.5mm	95mm	1
3-510	6.5mm	100mm	1
3-512	6.5mm	105mm	1
3-514	6.5mm	110mm	1
3-516	6.5mm	115mm	1
3-518	6.5mm	120mm	1

**Cancellous Screw 6.5mm (Stainless Steel)**  
Spongiosaschraube 6.5 (Stahl)

Screw Diameter • Schraubendurchmesser 6.50mm  
 Core Diameter • Kerndurchmesser 3.00mm  
 Head Diameter • Kopfdurchmesser 8.00mm  
 Drill Bit for threaded Hole • Bohrer für Kernloch 3.20mm



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-48000	6.5mm	25mm	1
3-48200	6.5mm	30mm	1
3-48400	6.5mm	35mm	1
3-48600	6.5mm	40mm	1
3-48800	6.5mm	45mm	1
3-49000	6.5mm	50mm	1
3-49200	6.5mm	55mm	1
3-49400	6.5mm	60mm	1
3-49600	6.5mm	65mm	1
3-49800	6.5mm	70mm	1
3-50000	6.5mm	75mm	1
3-50200	6.5mm	80mm	1
3-50400	6.5mm	85mm	1
3-50600	6.5mm	90mm	1
3-50800	6.5mm	95mm	1
3-51000	6.5mm	100mm	1
3-51200	6.5mm	105mm	1
3-51400	6.5mm	110mm	1
3-51600	6.5mm	115mm	1
3-51800	6.5mm	120mm	1

**Cancellous Screw 6.5mm (Titan)**  
Spongiosaschraube 6.5 (Titan)

Screw Diameter • Schraubendurchmesser 6.50mm  
 Core Diameter • Kerndurchmesser 3.00mm  
 Head Diameter • Kopfdurchmesser 8.00mm  
 Drill Bit for threaded Hole • Bohrer für Kernloch 3.20mm





Cat. no. Art. Nr.	∅	Material Material	PU VPE
3-522	13mm	Stainless Steel	1
3-5220	13mm	Titanium	1

Washer 13mm  
Unterlagscheibe 13mm



Cat. no. Art. Nr.	∅	Material Material	PU VPE
3-526	13mm	Stainless Steel	1
3-5260	13mm	Titanium	1

Washer with 4 Spikes 13mm  
Unterlagscheibe mit 4 Dornen 13mm



Cat. no. Art. Nr.	∅	Material Material	PU VPE
3-528	13mm	Stainless Steel	1
3-5280	13mm	Titanium	1

Washer with 10 Spikes 13mm  
Unterlagscheibe mit 10 Dornen 13mm



Cat. no. Art. Nr.	∅	Material Material	PU VPE
3-5201	11mm	Titanium	1
3-5220	13mm	Titanium	1

Washer (Titanium)  
Unterlagscheibe (Titan)



Cat. no. Art. Nr.	∅	Material Material	PU VPE
3-524	8mm	Stainless Steel	1

Nut for Cortical and Malleolar Screws 8mm (Stainless Steel)  
Mutter für Corticalis und Malleolarschrauben 8mm (Stahl)






Cat. no. Art. Nr.	∅	Length Länge	PU VPE
T 321	-	-	1

Aluminium Container (without Racks)  
Aluminium Container (ohne Racks)




Cat. no. Art. Nr.	for Screws ∅ für Schrauben ∅	PU VPE
T 310	1.5mm	1
T 311	2.0mm	1
T 312	2.7mm	1
T 313	3.5mm / 1.25 pitch	1
T 314	3.5mm / 1.75 pitch	1
T 315	4.0mm / half thread	1
T 316	4.0mm / full thread	1
T 317	4.5mm	1

Rack for Screws  
Rack für Schrauben



Cat. no. Art. Nr.	∅	Length Länge	PU VPE
3-530	-	100mm / 95mm	1


Measuring Gauge for Screws  
Messlehre für Schrauben





Set Cat. No. T 330

Basic Plate Set (Stainless Steel)  
Standart Plattenset (Stahl)



Listing: T 330 Plate Set Broad and Narrow Plates (Steel Plates an mat. DIN ISO 5832-1) for use dia 4.5mm and dia 6.5mm Screws



Set Cat. No. T 330

Cat. no. Art. Nr.	included items	QTY St.
T 331	Tray for Broad and Narrow Plates	1
<b>Semi-Tubular Plates</b>		
3-674	Semi-Tubular Plate 4 Holes 71mm	2
3-676	Semi-Tubular Plate 5 Holes 87mm	2
3-678	Semi-Tubular Plate 6 Holes 103mm	2
3-680	Semi-Tubular Plate 7 Holes 119mm	1
<b>Narrow dynamic compression Plates</b>		
3-604	Narrow DCP Plate 4 Holes 71mm	2
3-606	Narrow DCP Plate 5 Holes 87mm	2
3-608	Narrow DCP Plate 6 Holes 103mm	4
3-610	Narrow DCP Plate 7 Holes 119mm	2
3-612	Narrow DCP Plate 8 Holes 135mm	2
3-614	Narrow DCP Plate 9 Holes 151mm	2
3-616	Narrow DCP Plate 10 Holes 167mm	1
3-620	Narrow DCP Plate 12 Holes 199mm	1
<b>Broad dynamic compression Plates</b>		
3-640	Broad DCP Plate 6 Holes 103mm	2
3-642	Broad DCP Plate 7 Holes 119mm	2
3-644	Broad DCP Plate 8 Holes 135mm	2
3-646	Broad DCP Plate 9 Holes 151mm	1
3-648	Broad DCP Plate 10 Holes 167mm	1
3-652	Broad DCP Plate 12 Holes 199mm	1
3-656	Broad DCP Plate 14 Holes 231mm	1
<b>Spoon Plates</b>		
3-692	Spoon-Plate 5 Holes 100mm long	1
3-694	Spoon-Plate 6 Holes 120mm long	1
<b>T-Plates</b>		
3-702	T-Plate 4 Holes 84mm	1
3-706	T-Plate 6 Holes 116mm	1
3-708	T-Plate 8 Holes 148mm	1
<b>T-Buttress Plates</b>		
3-714	T-Buttress Plate 4 Holes 81mm	1
<b>L-Buttress Plates</b>		
3-720	L-Buttress Plate 4 Holes left	1
3-722	L-Buttress Plate 4 Holes right	1



Sterilization Container (recommended)

15-1870 Bottom non-perforated  
580 x 280 x 100 mm

15-1850 Lid perforated  
580 x 280mm silver



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-570	6	106mm	1
3-572	7	124mm	1
3-574	8	142mm	1
3-576	9	160mm	1
3-578	10	178mm	1
3-580	11	196mm	1
3-582	12	214mm	1
3-584	13	232mm	1
3-586	14	250mm	1
3-588	15	268mm	1
3-590	16	286mm	1
3-592	17	304mm	1
3-594	18	322mm	1

Broad MTP Plate (Stainless Steel)  
Breite MTP Platte (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-5700	6	106mm	1
3-5720	7	124mm	1
3-5740	8	142mm	1
3-5760	9	160mm	1
3-5780	10	178mm	1
3-5800	11	196mm	1
3-5820	12	214mm	1
3-5840	13	232mm	1
3-5860	14	250mm	1
3-5880	15	268mm	1
3-5900	16	286mm	1
3-5920	17	304mm	1
3-5940	18	322mm	1

Broad MTP Plate (Titanium)  
Breite MTP Platte (Titan)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-540	2	34mm	1
3-542	3	52mm	1
3-544	4	70mm	1
3-546	5	88mm	1
3-548	6	106mm	1
3-550	7	124mm	1
3-552	8	142mm	1
3-554	9	160mm	1
3-556	10	178mm	1
3-558	11	196mm	1
3-560	12	214mm	1
3-564	14	250mm	1
3-568	16	286mm	1

Narrow MTP Plate (Stainless Steel)  
Schmale MTP Platte (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-5400	2	34mm	1
3-5420	3	52mm	1
3-5440	4	70mm	1
3-5460	5	88mm	1
3-5480	6	106mm	1
3-5500	7	124mm	1
3-5520	8	142mm	1
3-5540	9	160mm	1
3-5560	10	178mm	1
3-5580	11	196mm	1
3-5600	12	214mm	1
3-5640	14	250mm	1
3-5680	16	286mm	1

Narrow MTP Plate (Titanium)  
Schmale MTP Platte (Titan)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-600	2	39mm	1
3-602	3	55mm	1
3-604	4	71mm	1
3-606	5	87mm	1
3-608	6	103mm	1
3-610	7	119mm	1
3-612	8	135mm	1
3-614	9	151mm	1
3-616	10	167mm	1
3-618	11	183mm	1
3-620	12	199mm	1
3-622	13	215mm	1
3-624	14	231mm	1
3-626	15	247mm	1
3-628	16	263mm	1

Narrow DCP Plate (Stainless Steel)  
Schmale DCP Platte (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-640	6	103mm	1
3-642	7	119mm	1
3-644	8	135mm	1
3-646	9	151mm	1
3-648	10	167mm	1
3-650	11	183mm	1
3-652	12	199mm	1
3-654	13	215mm	1
3-656	14	231mm	1
3-658	15	247mm	1
3-660	16	263mm	1
3-662	17	279mm	1
3-664	18	295mm	1

Broad DCP Plate (Stainless Steel)  
Breite DCP Platte (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-670	2	39mm	1
3-672	3	55mm	1
3-674	4	71mm	1
3-676	5	87mm	1
3-678	6	103mm	1
3-680	7	119mm	1
3-682	8	135mm	1
3-684	9	151mm	1
3-686	10	167mm	1
3-688	11	183mm	1
3-690	12	199mm	1

Sami-Tubular Plate (Stainless Steel)  
Halbrohrplatte (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-6700	2	39mm	1
3-6720	3	55mm	1
3-6740	4	71mm	1
3-6760	5	87mm	1
3-6780	6	103mm	1
3-6800	7	119mm	1
3-6820	8	135mm	1
3-6840	9	151mm	1
3-6860	10	167mm	1
3-6880	11	183mm	1
3-6900	12	199mm	1

Sami-Tubular Plate (Titanium)  
Halbrohrplatte (Titan)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-692	5	100mm	1
3-694	6	120mm	1

Spoon Plate (Stainless Steel)  
Löffelplatte (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-6920	5	100mm	1
3-6940	6	120mm	1

Spoon Plate (Titanium)  
Löffelplatte (Titan)





Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-700	3	68mm	1
3-702	4	84mm	1
3-704	5	100mm	1
3-706	6	116mm	1
3-708	8	148mm	1

T-Plate (Stainless Steel)  
T-Platte (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-7000	3	68mm	1
3-7020	4	84mm	1
3-7040	5	100mm	1
3-7060	6	116mm	1
3-7080	8	148mm	1

T-Plate (Titanium)  
T-Platte (Titan)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-714	4	81mm	1
3-716	5	96mm	1
3-718	6	112mm	1

T-Buttress Plate (Stainless Steel)  
T-Abstützplatte (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-7140	4	81mm	1
3-7160	5	96mm	1
3-7180	6	112mm	1

T-Buttress Plate (Titanium)  
T-Abstützplatte (Titan)




Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-720	4 left	-	1
3-722	4 right	-	1
3-724	5 left	-	1
3-726	5 right	-	1
3-728	6 left	-	1
3-730	6 right	-	1

L-Buttress Plate for Tibia (Stainless Steel)  
L-Abstützplatte für Tibia (Stahl)




Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-7200	4 left	-	1
3-7220	4 right	-	1
3-7240	5 left	-	1
3-7260	5 right	-	1
3-7280	6 left	-	1
3-7300	6 right	-	1

L-Buttress Plate (Titanium)  
L-Abstützplatte (Titan)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-750	8	135mm	1
3-752	8	145mm	1
3-754	8	155mm	1
3-756	8	165mm	1

Broad Lengthening Plate 8 Holes (Stainless Steel)  
Breite Verlängerungsplatte 8 Loch (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-760	8	135mm	1
3-762	8	145mm	1
3-764	8	155mm	1
3-766	8	165mm	1
3-768	8	175mm	1
3-770	8	185mm	1

Narrow Lengthening Plate 8 Holes (Stainless Steel)  
Schmale Verlängerungsplatte 8 Loch (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-772	10	179mm	1
3-774	10	189mm	1
3-776	10	199mm	1
3-778	10	209mm	1
3-779	10	219mm	1
3-780	10	229mm	1
3-782	10	239mm	1
3-784	10	249mm	1

Broad Lengthening Plate 10 Holes (Stainless Steel)  
Breite Verlängerungsplatte 10 Loch (Stahl)




Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-820	8	170mm	1
3-822	9	186mm	1
3-824	10	202mm	1
3-826	11	218mm	1

DCP Cobra Plate (Stainless Steel)  
DCP Kreuzplatte (Stahl)




Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-800	7 left	158mm	1
3-802	7 right	158mm	1
3-804	9 left	190mm	1
3-806	9 right	190mm	1
3-808	11 left	221mm	1
3-810	11 right	221mm	1
3-812	13 left	253mm	1
3-814	13 right	253mm	1
3-816	15 left	285mm	1
3-818	15 right	285mm	1

Condylar Buttress Plate left and right (Stainless Steel)  
Kondylenabstützplatte links und rechts (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-8000	7 left	158mm	1
3-8020	7 right	158mm	1
3-8040	9 left	190mm	1
3-8060	9 right	190mm	1
3-8080	11 left	221mm	1
3-8100	11 right	221mm	1
3-8120	13 left	253mm	1
3-8140	13 right	253mm	1
3-8160	15 left	285mm	1
3-8180	15 right	285mm	1


Condylar Buttress Plate left and right (Titanium)  
Kondylenabstützplatte links und rechts (Titan)






Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-830	5 right	-	1
3-831	5 left	-	1
3-832	7 right	-	1
3-833	7 left	-	1
3-834	9 right	-	1
3-835	9 left	-	1
3-836	11 right	-	1
3-837	11 left	-	1
3-838	13 right	-	1
3-839	13 left	-	1

Lateral Tibia Plate left and right (Stainless Steel)  
Laterale Tibiaplatte links und rechts (Stahl)




Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-8300	5 right	-	1
3-8310	5 left	-	1
3-8320	7 right	-	1
3-8330	7 left	-	1
3-8340	9 right	-	1
3-8350	9 left	-	1
3-8360	11 right	-	1
3-8370	11 left	-	1
3-8380	13 right	-	1
3-8390	13 left	-	1

Lateral Tibia Plate left and right (Titanium)  
Laterale Tibiaplatte links und rechts (Titan)




Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-850	7 right	122mm	1
3-852	7 left	122mm	1
3-854	9 right	158mm	1
3-856	9 left	158mm	1
3-858	11 right	194mm	1
3-860	11 left	194mm	1
3-862	13 right	230mm	1
3-866	13 left	230mm	1

Distal Tibia Plate left and right (Stainless Steel)  
Distale Tibiaplatte links und rechts (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-8500	7 right	122mm	1
3-8520	7 left	122mm	1
3-8540	9 right	158mm	1
3-8560	9 left	158mm	1
3-8580	11 right	194mm	1
3-8600	11 left	194mm	1
3-8620	13 right	230mm	1
3-8660	13 left	230mm	1

Distal Tibia Plate left and right (Titanium)  
Distale Tibiaplatte links und rechts (Titan)




Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-870	7 right	122mm	1
3-872	7 left	122mm	1
3-874	9 right	158mm	1
3-876	9 left	158mm	1
3-878	11 right	194mm	1
3-880	11 left	194mm	1
3-882	13 right	230mm	1
3-884	13 left	230mm	1

Proximal Tibia Plate left and right (Stainless Steel)  
Proximale Tibiaplatte links und rechts (Stahl)




Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-8700	7 right	122mm	1
3-8720	7 left	122mm	1
3-8740	9 right	158mm	1
3-8760	9 left	158mm	1
3-8780	11 right	194mm	1
3-8800	11 left	194mm	1
3-8820	13 right	230mm	1
3-8840	13 left	230mm	1

Proximal Tibia Plate left and right (Titanium)  
Proximale Tibiaplatte links und rechts (Titan)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-886	7 right	122mm	1
3-888	7 left	122mm	1
3-890	9 right	158mm	1
3-892	9 left	158mm	1
3-894	11 right	194mm	1
3-896	11 left	194mm	1
3-898	13 right	230mm	1
3-900	13 left	230mm	1

Distal Femur Plate left and right (Stainless Steel)  
Distale Femurplatte links und rechts (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-8860	7 right	122mm	1
3-8880	7 left	122mm	1
3-8900	9 right	158mm	1
3-8920	9 left	158mm	1
3-8940	11 right	194mm	1
3-8960	11 left	194mm	1
3-8980	13 right	230mm	1
3-9000	13 left	230mm	1


Distal Femur Plate left and right (Titanium)  
Distale Femurplatte links und rechts (Titan)






Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-902	7 right	122mm	1
3-904	7 left	122mm	1
3-906	9 right	158mm	1
3-908	9 left	158mm	1
3-910	11 right	194mm	1
3-912	11 left	194mm	1
3-914	13 right	230mm	1
3-916	13 left	230mm	1

Distal Tibia Plate for Fibula left and right (Stainless Steel)  
Distale Tibiaplatte für Fibula links und rechts (Stahl)



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-9020	7 right	122mm	1
3-9040	7 left	122mm	1
3-9060	9 right	158mm	1
3-9080	9 left	158mm	1
3-9100	11 right	194mm	1
3-9120	11 left	194mm	1
3-9140	13 right	230mm	1
3-9160	13 left	230mm	1

Distal Tibia Plate for Fibula left and right (Titanium)  
Distale Tibiaplatte für Fibula links und rechts (Titan)







TREU-Instrumente GmbH  
take-off Gewerbepark 130-132  
78579 Neuhausen ob Eck  
Germany

Phone +49 (0) 74 67/947 673 0  
Fax +49 (0) 74 67/947 673 50

[www.treu.com](http://www.treu.com)  
[info@treu.com](mailto:info@treu.com)