






YOUR RELIABLE PARTNER IN TRAUMA CARE

WE FOCUS THE FUTURE

Pelvic Instruments and Implants







Set Cat. No. T 1700


Cat. no. Art. Nr.	included items	QTY St.
15-1870	Bottom non-perforated 580x280x100mm	1
15-1850	Lid perforated 580x280mm silver	1
15-3012	Wire Basket Stainless Steel 540x255x50mm	1
9-021	Drill Bit AO Shaft Ø 3.2mm x 195/170mm	2
9-027	Drill Bit AO Shaft Ø 4.5mm x 195/170mm	2
9-074	Tap Thread Pitch Ø 3.5mm dia 1.25 AO	1
9-078	Tap AO Shaft long Thread Ø 4.5mm	1
9-142	Hexagonal Screwdriver 2.5mm	1
9-144	Hexagonal Screwdriver	1
9-114	Depth Gauge for Large Screws	1
17-1400	Bending Pliers for 2.7mm/3.5mm Plates	1
17-1424	Bending Templates 3.5mm curved	1
17-1426	Bending Templates 3.5mm straight	1
17-1404	Reposition Forceps	1

Palvic Basic Set I
Beckenrepositionsinstrumentarium I



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
17-1400	-	265mm	1

Bending Pliers for 2.7mm/3.5mm Plates
Biegezange für 2.7mm/3.5mm Platten




Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
17-1402	-	195mm	1

Bending Iron for 3.5mm/4.5mm Plates
Biegeeisen für 3.5mm/4.5mm Platten



Cat. no. Art. Nr.	Ø		PU VPE
17-1420	4.5mm	curved	1
17-1422	4.5mm	straight	1
17-1424	3.5mm	curved	1
17-1426	3.5mm	straight	1

Bending Templates
Biegeschablone



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
17-1430	15.0mm	220mm	1
17-1432	20.0mm	220mm	1

Osteotomy Chisel
Osteotomie Meissel





Set Cat. No. T 1710


Cat. no. Art. Nr.	included items	QTY St.
15-1870	Bottom non-perforated 580x280x100mm	1
15-1850	Lid perforated 580x280mm silver	1
15-3014	Wire Basket Stainless Steel 540x255x70mm	1
15-3266-2-1	Standard Silicone Device	4
9-090	T-Handle	1
11-4604	Sharp Bone Hook 230mm	1
9-676	Retractor widetip 24.0mm	1
9-521	Pelvic Reposition Forceps	1
9-522	Farabeuf-Lambotte Forceps 260mm	1
9-604	Reduction Forceps 200mm	1
17-1550	Pelvic Reposition Forceps 190mm	1
17-1554	Pelvic Reposition Forceps 230mm	1
17-1552	Pelvic Reposition Forceps 250mm	1
17-1556	Pelvic Reposition Forceps 400mm	1
17-1558	Pelvic Reposition Forceps 400mm	1
17-1600	Straight Ball Spike	1
8-1516	Schanz Screw 6.0x200mm	4

Palvic Basic Set II
Beckenrepositionsinstrumentarium II



Cat. no. Art. Nr.	∅	Length Länge	PU VPE
9-090	-	130mm	1

T-Handle (9-089) with Universal three-Jaw Chuck (9-091) and Key (9-0910) com.
T-Griff (9-089) mit Universalfutter (9-091) und Schlüssel (9-0910) komplett



Cat. no. Art. Nr.	∅	Length Länge	PU VPE
11-4604	-	230mm	1

Sharp Bone Hook 230mm
Einzinkhaken mit Quergriff 230mm



Cat. no. Art. Nr.	∅	Length Länge	PU VPE
9-676	-	270mm	1


Retractor widetip 24.0mm
Knochenheber breite Spitze 24.0mm





Cat. no. Art. Nr.	∅	Length Länge	PU VPE
9-521	3.5mm	190mm	1
9-523	4.5mm	190mm	1

Pelvic Reposition Forceps
Becken Repositionszange




Cat. no. Art. Nr.	∅	Length Länge	PU VPE
9-604	-	200mm	1

Reduction Forceps 200mm
Reduktionszange 200mm




Cat. no. Art. Nr.	∅	Length Länge	PU VPE
17-1550	-	200mm	1

Pelvic Reposition Forceps 200mm
Becken Repositionszange 200mm




Cat. no. Art. Nr.	∅	Length Länge	PU VPE
17-1554	-	245mm	1

Pelvic Reposition Forceps 245mm
Becken Repositionszange 245mm



Cat. no. Art. Nr.	∅	Length Länge	PU VPE
17-1552	-	255mm	1

Pelvic Reposition Forceps 255mm
Becken Repositionszange 255mm



Cat. no. Art. Nr.	∅	Length Länge	PU VPE
17-1556	-	400mm	1

Pelvic Reposition Forceps 400mm
Becken Repositionszange 400mm



Cat. no. Art. Nr.	∅	Length Länge	PU VPE
17-1558	-	405mm	1

Pelvic Reposition Forceps 405mm
Becken Repositionszange 405mm



Cat. no. Art. Nr.	∅	Length Länge	PU VPE
17-1600	-	300mm	1

Straight Ball Spike
Gerade Kugelspitze




Cat. no. Art. Nr.	∅	Length Länge	PU VPE
8-1516	6.0mm	200mm	1

Schanz Screw 6.0x200mm
Schanzschraube 6.0x200mm



Cat. no. Art. Nr.	∅	Length Länge	PU VPE
17-1404	-	340mm	1

Reposition Forceps
Repositionszange





Set Cat. No. T 1720

Cat. no. Art. Nr.	included items	QTY St.
T 111	Aluminium Case	1
T 112	Tray	1
T 402	Lower Tray	1
3-600	Narrow DCP Plate 2 Holes 39mm	2
17-1100	Reconstruction Plate 4.5/45mm 3 Holes	2
17-1102	Reconstruction Plate 4.5/61mm 4 Holes	2
17-1122	Reconstruction Plate 4.5/221mm 14 Holes	1
4-270	Reconstruction Plate 3.5mm 5 Holes	2
4-271	Reconstruction Plate 3.5mm 6 Holes	2
4-272	Reconstruction Plate 3.5mm 7 Holes	2
4-273	Reconstruction Plate 3.5mm 8 Holes	2
4-274	Reconstruction Plate 3.5mm 9 Holes	2
4-275	Reconstruction Plate 3.5mm 10 Holes	1
4-277	Reconstruction Plate 3.5mm 12 Holes	1
4-279	Reconstruction Plate 3.5mm 14 Holes	1
4-281	Reconstruction Plate 3.5mm 16 Holes	1
4-282	Reconstruction Plate 3.5mm 18 Holes	1
17-1140	Reconstruction Plate 3.5/70mm 6 Holes	2
17-1142	Reconstruction Plate 3.5/94mm 8 Holes	2
17-1144	Reconstruction Plate 3.5/118mm 10 Holes	1
17-1146	Reconstruction Plate 3.5/142mm 12 Holes	1
17-1148	Reconstruction Plate 3.5/166mm 14 Holes	1
17-1150	Reconstruction Plate 3.5/190mm 16 Holes	1
3-1250	Cortical Screw 4.5 x 70mm	3
3-1252	Cortical Screw 4.5 x 75mm	3
3-1254	Cortical Screw 4.5 x 80mm	3
3-1256	Cortical Screw 4.5 x 85mm	3
3-1258	Cortical Screw 4.5 x 90mm	3
3-1260	Cortical Screw 4.5 x 95mm	3
3-1262	Cortical Screw 4.5 x 100mm	3
3-1264	Cortical Screw 4.5 x 105mm	3
3-1266	Cortical Screw 4.5 x 110mm	3
4-1034	Cortical Screw 3.5 x 55mm	3
4-1036	Cortical Screw 3.5 x 60mm	3
4-1038	Cortical Screw 3.5 x 65mm	3
4-1040	Cortical Screw 3.5 x 70mm	3
4-1042	Cortical Screw 3.5 x 75mm	3
4-1044	Cortical Screw 3.5 x 80mm	3
3-415	Spongiosa Screw 6.5 x 115mm	3
3-416	Spongiosa Screw 6.5 x 120mm	3
3-490	Spongiosa Screw 6.5 x 50mm full Thread	3
3-492	Spongiosa Screw 6.5 x 55mm full Thread	3
3-494	Spongiosa Screw 6.5 x 60mm full Thread	3

Palvic Basic Implant Set
Becken-Repositionsimplantatesatz





Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
17-1100	3	45mm	1
17-1102	4	61mm	1
17-1104	5	77mm	1
17-1106	6	93mm	1
17-1108	7	109mm	1
17-1110	8	125mm	1
17-1112	9	141mm	1
17-1114	10	157mm	1
17-1116	11	173mm	1
17-1118	12	189mm	1
17-1120	13	205mm	1
17-1122	14	221mm	1
17-1124	15	237mm	1
17-1126	16	253mm	1

Reconstruction Plate 4.5mm (Stainless Steel)
Rekonstruktionsplatte 4.5mm (Stahl)

Profil • Profil 12.0 x 2.8mm
Distance between Holes • Abstand zwischen den Löchern 16.0mm



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
17-11000	3	45mm	1
17-11020	4	61mm	1
17-11040	5	77mm	1
17-11060	6	93mm	1
17-11080	7	109mm	1
17-11100	8	125mm	1
17-11120	9	141mm	1
17-11140	10	157mm	1
17-11160	11	173mm	1
17-11180	12	189mm	1
17-11200	13	205mm	1
17-11220	14	221mm	1
17-11240	15	237mm	1
17-11260	16	253mm	1

Reconstruction Plate 4.5mm (Titanium)
Rekonstruktionsplatte 4.5mm (Titan)

Profil • Profil 12.0 x 2.8mm
Distance between Holes • Abstand zwischen den Löchern 16.0mm





Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
4-270	5	58mm	1
4-271	6	70mm	1
4-272	7	82mm	1
4-273	8	94mm	1
4-274	9	106mm	1
4-275	10	118mm	1
4-276	11	130mm	1
4-277	12	142mm	1
4-278	13	154mm	1
4-279	14	166mm	1
4-280	15	178mm	1
4-281	16	190mm	1
4-282	18	214mm	1
4-283	20	238mm	1
4-284	22	262mm	1

Reconstruction Plate 3.5mm (Stainless Steel)
Rekonstruktionsplatte 3.5mm (Stahl)

Profil · Profil

12.0 x 2.8mm

Distance between Holes · Abstand zwischen den Löchern

16.0mm



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
4-2700	5	58mm	1
4-2710	6	70mm	1
4-2720	7	82mm	1
4-2730	8	94mm	1
4-2740	9	106mm	1
4-2750	10	118mm	1
4-2760	11	130mm	1
4-2770	12	142mm	1
4-2780	13	154mm	1
4-2790	14	166mm	1
4-2800	15	178mm	1
4-2810	16	190mm	1
4-2820	18	214mm	1
4-2830	20	238mm	1
4-2840	22	262mm	1

Reconstruction Plate 3.5mm (Titanium)
Rekonstruktionsplatte 3.5mm (Titan)

Profil · Profil

12.0 x 2.8mm

Distance between Holes · Abstand zwischen den Löchern

16.0mm





Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
17-1140	6	70mm	1
17-1142	8	94mm	1
17-1144	10	118mm	1
17-1146	12	142mm	1
17-1148	14	166mm	1
17-1150	16	190mm	1
17-1152	18	214mm	1

Curved Reconstruction Plate 3.5mm (Stainless Steel)
Gebogene Rekonstruktionsplatte 3.5mm (Stahl)

Profil • Profil 10.0 x 2.8mm
 Distance between Holes • Abstand zwischen den Löchern 12.0mm
 Curving Radius • Krümmungsradius 100.0mm



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
17-11400	6	70mm	1
17-11420	8	94mm	1
17-11440	10	118mm	1
17-11460	12	142mm	1
17-11480	14	166mm	1
17-11500	16	190mm	1
17-11520	18	214mm	1

Curved Reconstruction Plate 3.5mm (Titanium)
Gebogene Rekonstruktionsplatte 3.5mm (Titan)

Profil • Profil 10.0 x 2.8mm
 Distance between Holes • Abstand zwischen den Löchern 12.0mm
 Curving Radius • Krümmungsradius 100.0mm



Cat. no. Art. Nr.	Holes Löcher	Length Länge	PU VPE
3-600	2	39mm	1

Narrow DCP Plate 2 Holes, 39mm length
Schmale DCP-Platte 2 Loch, 39mm Länge






Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-1200	4.5mm	14mm	1
3-1202	4.5mm	16mm	1
3-1204	4.5mm	18mm	1
3-1206	4.5mm	20mm	1
3-1208	4.5mm	22mm	1
3-1210	4.5mm	24mm	1
3-1212	4.5mm	26mm	1
3-1214	4.5mm	28mm	1
3-1216	4.5mm	30mm	1
3-1218	4.5mm	32mm	1
3-1220	4.5mm	34mm	1
3-1222	4.5mm	36mm	1
3-1224	4.5mm	38mm	1
3-1226	4.5mm	40mm	1
3-1228	4.5mm	42mm	1
3-1230	4.5mm	44mm	1
3-1232	4.5mm	46mm	1
3-1234	4.5mm	48mm	1
3-1236	4.5mm	50mm	1
3-1238	4.5mm	52mm	1
3-1240	4.5mm	54mm	1
3-1242	4.5mm	56mm	1
3-1244	4.5mm	58mm	1
3-1246	4.5mm	60mm	1
3-1248	4.5mm	65mm	1
3-1250	4.5mm	70mm	1
3-1252	4.5mm	75mm	1
3-1254	4.5mm	80mm	1
3-1256	4.5mm	85mm	1
3-1258	4.5mm	90mm	1
3-1260	4.5mm	95mm	1
3-1262	4.5mm	100mm	1
3-1264	4.5mm	105mm	1
3-1266	4.5mm	110mm	1
3-1268	4.5mm	115mm	1
3-1270	4.5mm	120mm	1

Cortical Screw 4.5mm, self-cutting (Stainless Steel)
Kortikalisschraube 4.5mm, selbstschneidend (Stahl)

Screw Diameter • Schraubendurchmesser 4.50mm
Core Diameter • Kerndurchmesser 3.00mm
Head Diameter • Kopfdurchmesser 8.00mm



self-cutting Screws, tapping not needed



Cat. no. Art. Nr.	Ø	Length Länge	PU VPE
3-1280	4.5mm	22mm	1
3-1282	4.5mm	24mm	1
3-1284	4.5mm	26mm	1
3-1286	4.5mm	28mm	1
3-1288	4.5mm	30mm	1
3-1290	4.5mm	32mm	1
3-1292	4.5mm	34mm	1
3-1294	4.5mm	36mm	1
3-1296	4.5mm	38mm	1
3-1298	4.5mm	40mm	1
3-1300	4.5mm	42mm	1
3-1302	4.5mm	44mm	1
3-1304	4.5mm	46mm	1
3-1306	4.5mm	48mm	1
3-1308	4.5mm	50mm	1
3-1310	4.5mm	52mm	1
3-1312	4.5mm	54mm	1
3-1314	4.5mm	56mm	1
3-1316	4.5mm	58mm	1
3-1318	4.5mm	60mm	1
3-1320	4.5mm	65mm	1
3-1322	4.5mm	70mm	1
3-1324	4.5mm	75mm	1
3-1326	4.5mm	80mm	1
3-1328	4.5mm	85mm	1
3-1330	4.5mm	90mm	1
3-1332	4.5mm	95mm	1
3-1334	4.5mm	100mm	1
3-1336	4.5mm	105mm	1
3-1338	4.5mm	110mm	1

Shank Screw 4.5mm, self-cutting (Stainless Steel)
 Schaftschraube 4.5mm, selbstschneidend (Stahl)

Screw Diameter • Schraubendurchmesser 4.50mm
 Shaft Diameter • Schaftdurchmesser 4.50mm
 Head Diameter • Kopfdurchmesser 8.00mm



self-cutting Screws, tapping not needed





TREU-Instrumente GmbH
take-off Gewerbepark 130-132
78579 Neuhausen ob Eck
Germany

Phone +49 (0) 74 67/947 673 0
Fax +49 (0) 74 67/947 673 50

www.treu.com
info@treu.com