

Canula / Biopsy needle with stylet

Ref.-no. EMS 23-10
Working length 100 mm
Diameter Ø 0.65 mm / 23 G



- Sharp, atraumatic and smooth
- Fixation stylet to prevent unwanted puncture infection
- Special cutting point of needle tip and specialized surface treatment allow easy puncturing even of deeper tissue
- For better positioning the needles are marked with a galvanically applied unremovable depth gauge at intervals of 1 cm
- Special stainless steel stylets provide precise access to deeper tissue regions. This material guarantees extreme rigidity
- Special, flexible stainless steel for easier positioning correction, granting high stability

Micromed Medizintechnik GmbH
Eisenbahnstr. 84
D-78573 Wurmlingen

Tel +49 7461 9648 550
Fax +49 7461 9648 5595

info@micromed.com
www.micromed.com