

## MicroPen<sup>®</sup>

for RF-Units of  
Micromed  
ERBE  
Martin  
Berchtold  
Aesculap  
Valleylab



### Cable

- in various lengths 3m, 4m, 5m
- with various types of plugs for different manufacturers' RF surgical units

### Handle and cable separable

- flexible handle designed to fit into different RF-generators
- interchangeable cables

### Electrode handle MicroPen


- aesthetic design
- optimum handling
- long life

### Push button

- lift of stroke can be felt easily
- click control
- high security by durable push button
- push button completely covered with insulating silicone

### Electrodes


- hexagonal fit for electrodes measuring 4 mm and 3/32 inches in diameter
- circular fitting for 1,6 mm microelectrodes
- electrodes can be exchanged easily



You want to know more?

To get further information  
please visit [www.micromed.biz](http://www.micromed.biz)

  
MICROMED®



100-104-016

100-201-040

100-201-024

100-101-040

100-203-040

MICROMED MICROPEN

**12**

12 MONTH GUARANTEE\*

#### Optimum synthetic substance

- high voltage strength up to 7500 Volt
- great number of sterilization cycles
- highly resistant to detergents

Electrode handle MicroPen® • Highest quality with 12 month guarantee • Electrode handle MicroPen® from Micromed®

The **MicroPen®** is an innovative electrode handle that can live up to highest demands. Due to optimum technical production and the use of top quality materials we have created a handle with maximum longevity as well as optimum handling comforts.

#### Handle and cable separable

By separating handle and cable the components can be replaced separately in case of defects. This lowers costs and is environmentally friendly.

#### Highly convenient

Sensibly located pushbuttons and the easily felt pressure point (on/off click) as well as the rough outside surface make the MicroPen especially user-friendly.

#### Cleaning and Disinfection

The MicroPen can be sterilized more than 250 times.


#### Highest compatibility

The MicroPen is compatible to a wide range of different manufacturers' surgical units, as for instance: Aesculap, Berchtold, Integra, Erbe, Martin and Micromed.

\*every MicroPen® comes with a 1 year guarantee on purchase.



[www.micromed.biz](http://www.micromed.biz)



Micromed Medizintechnik GmbH  
Esslinger Str. 37  
D-78532 Tuttlingen

Tel +49 7462 923250  
Fax +49 7462 923246  
[info@micromed.biz](mailto:info@micromed.biz)